

## Drug Analysis

### 1 Scope

This procedure is used for the analysis of items suspected of containing drugs. Questioned items may consist of a variety of substances including unknown solids and liquids, prescription medications, and over-the-counter products [see *Analysis of Tablets and Capsules (GenChem 7)* for prescription and over-the-counter tablets, capsules, etc.]. While the below procedure describes many techniques, identification of a drug relies upon positive results from two orthogonal techniques with at least one of the techniques providing structural elucidation information.

This procedure applies to Chemistry Unit (CU) personnel that are qualified and authorized to examine General Chemistry evidence for the presence of drugs.

### 2 Equipment/Materials/Reagents

- Common laboratory glassware and equipment
- Analytical balance
- Digital microscope
- Stereo microscope
- Ultraviolet light source (long wavelength)
- Acetaldehyde
- Acetonitrile
- Chloroform
- Cobalt thiocyanate
- Deionized water
- Diethyl ether
- Formaldehyde (40%)
- Hydrochloric acid
- Methanol (MeOH)
- Nitric acid
- Sodium bicarbonate
- Sodium carbonate
- Sodium nitroprusside (aka sodium nitroferricyanide)
- Sodium sulfate (anhydrous)
- Sulfuric acid

- Evaporator
- Fourier Transform Infrared (FTIR) spectrophotometer with Attenuated Total Reflectance (ATR) or microscope attachment
- Polyethylene glycol (PEG, 550 average molecular weight)
- Time-of-flight mass spectrometer with direct analysis in real time ionization source (DART/TOFMS)
- Gas chromatograph/mass spectrometer (GC/MS) equipped with electron impact ionization and a 30 meter DB-5 column (or equivalent)
- Gas chromatograph/mass spectrometer (GC/MS) equipped with chemical ionization and a 30 meter DB-5 column (or equivalent)

### **3 Standards and Controls**

#### **3.1 Negative Control**

The same volume of solvent from the same source and lot used to extract or rinse the questioned item(s) and within a similar container (e.g., test tube, vial) will be used as the Negative Control.

#### **3.2 Positive Control**

Prepared by making a 1 mg/mL (as base) stock standard solution of the drug within a suitable solvent. A working standard solution of 0.1 mg/mL is typically used with GC/MS by diluting the stock standard 1:10. These solutions will be stored in a freezer or refrigerator. Other concentrations may be prepared and used as needed.

### **4 Preparation of Color Test Reagents**

#### **4.1 Marquis Reagent**

Prepared by adding 8-10 drops of 40% formaldehyde to 10 mL of concentrated sulfuric acid. This solution is stored at room temperature in an amber glass bottle. Discard the solution when it begins to discolor.

## 4.2 Scott's Reagent

- Reagent A- Prepared by adding 2 grams of cobalt thiocyanate to 100 mL of deionized water. Mix thoroughly. The solution should be pink in color. This solution is stored in a glass bottle at room temperature.
- Reagent B- Prepared by adding 8.5 mL of concentrated hydrochloric acid to 80 mL of deionized water and then diluting to 100 mL with deionized water. This solution is stored in a glass bottle at room temperature.

## 4.3 Sodium Nitroprusside Reagent

- Reagent A- Prepared by dissolving 1.1 grams of sodium nitroprusside (aka sodium nitroferricyanide) into 100 mL of deionized water and 4 mL of acetaldehyde. This solution is stored in an amber glass bottle in a refrigerator.
- Reagent B- Prepared fresh by dissolving 2 grams of sodium carbonate in 100 mL of deionized water. A small amount of solid sodium carbonate (or sodium bicarbonate) can be used in lieu of the aqueous solution.

## 5 Sampling

Sampling is performed according to the *Sampling Guidelines for Bulk Materials and Multi-Unit Populations*– General Chemistry SOP manual.

When non-statistical sampling is utilized on a heterogeneous item, the results of examinations will be clearly limited to the sample(s) that were selected and examined.

## 6 Procedure

Refer to *General Chemistry Instrument Parameters* (GenChem 34) for specific instrument settings and decision criteria.

- a. Use a traceable analytical balance to record the weight for each item, as relevant. For example, the weight of a drug paraphernalia item prior to rinsing is not typically relevant.
- b. Perform a visual examination of each item. Microscopy may be used as deemed necessary. Suspected lysergic acid diethylamide (LSD) samples may be observed under long wavelength ultraviolet light for blue fluorescence.

- c. Color tests may be performed. All samples (e.g., controls, items) will be added to the reagent to ensure that the spot plate well is free of contamination. Test tubes or other containers may be used in place of spot plates. The reagent and sample amounts may be adjusted as necessary. Color tests not described below may be prepared and used following similar practices. Include a copy of the reference relied upon for the color test. Examples of resources include *Clarke's Analysis of Drugs and Poisons*, SWGDRUG Drug Monographs, and the DEA *Analysis of Drugs Manual*.

#### Marquis Test:

- Add 2-3 drops of Marquis Reagent to the required number of spot plate wells. Add samples to separate wells and observe and record any changes. One well will remain unaltered during the exam to demonstrate the Marquis Reagent does not change color spontaneously. A purple color indicates the possible presence of an opiate. A violet to black color indicates the possible presence of 3,4-methylenedioxyamphetamine (MDA) or 3,4-methylenedioxymethamphetamine (MDMA). An orange color indicates the possible presence of an amphetamine compound.

#### Nitroprusside Test:

- Add 2-3 drops of Sodium Nitroprusside Reagent A to the required number of spot plate wells. Add samples to separate wells and observe and record any changes. Next, add 2-3 drops of Sodium Nitroprusside Reagent B (or a small amount of solid sodium carbonate or sodium bicarbonate) to each well and observe and record any changes. One well will only have Reagents A and B added to it to demonstrate that a color change does not occur. A blue or violet color upon addition of Reagent B indicates the possible presence of an amphetamine compound.

#### Scott's Test:

- Add 2-3 drops of Scott's Reagent A to the required number of spot plate wells. Add samples to separate wells and observe and record any changes. Next, add 2-3 drops of Scott's Reagent B to each well and observe and record any changes. One well will only have Reagents A and B added to it to demonstrate that a color change does not occur. A blue color prior to the addition of Reagent B indicates the possible presence of cocaine hydrochloride, whereas a blue color after the addition of Reagent B indicates the possible presence of cocaine base.

#### Concentrated Nitric Acid Test:

- Add 2-3 drops of concentrated nitric acid to the required number of spot plate wells. Add samples to separate wells and observe and record any changes. One well will remain unaltered during the exam to demonstrate the nitric acid does not change color spontaneously. An orange color that fumes indicates the possible presence of acetaminophen. An orange color that does not fume indicates the possible presence of morphine or codeine. A lime green color indicates the possible presence of guaifenesin or methocarbamol. A blue fluorescence under ultraviolet light indicates the possible presence of quinine.
- d. A representative sample may be analyzed by FTIR-ATR. A solid/powder item may be homogenized using a mortar and pestle. Liquid samples may be analyzed neat and/or allowed to evaporate on the ATR cell. The FTIR microscope attachment may be used as appropriate.
- e. If the item is a mixture, the following basic sequential solvent extraction may be used to isolate individual compounds for further analysis. Most basic organic drugs will be soluble in diethyl ether, while most drug salts will be soluble in chloroform. Sugars will usually be soluble in methanol. Utilize a fume hood for this process.
- Homogenize a representative portion of the item with a mortar and pestle.
  - Place a piece of folded filter paper into a funnel and place the funnel over an evaporating dish. Place the homogenized powder in the filter paper.
  - Pour 2 to 3 mL of diethyl ether over the powder and collect the filtrate into the evaporating dish. Set the dish to the side and allow the diethyl ether extract to evaporate.
  - Place a new evaporating dish under the funnel and wash any remaining powder with 2 to 3 mL of chloroform. Set the dish to the side and allow the chloroform extract to evaporate.
  - Place a new evaporating dish under the funnel and wash any remaining powder with 2 to 3 mL of methanol. Allow the methanol extract to evaporate.
  - Any remaining solids in the filter paper, as well as any solids recovered from the evaporated extracts may be analyzed by FTIR-ATR.
- f. An appropriate amount of the item may be placed into a small test tube and extracted or dissolved in an appropriate solvent (e.g., methanol, chloroform, acetonitrile) to achieve the desired concentration. For example, a concentration of ~1 mg/mL may be desirable for DART/TOFMS analysis, while a concentration of ~100 ug/mL is common for GC/MS analysis. If necessary, filter the solution or centrifuge and decant to remove any undissolved particulates. The Negative Control will be filtered or centrifuged as well. It may be necessary

to utilize acidic or basic conditions to more efficiently extract some drugs, see Appendix A for acid/neutral and alkaline drug extraction steps.

- g. If the item is an empty syringe (or similar item), rinse the interior of the barrel with a volume of an appropriate solvent (e.g., methanol, chloroform, acetonitrile) that is approximately equivalent to the volume capacity of the barrel. If the needle is intact and needs to be preserved for future DNA exams, then remove the plunger (if present) to introduce the rinse solvent to the barrel. Use the rinse to sample the barrel multiple times, then transfer the rinse to a labeled test tube. If necessary, concentrate the rinse and Negative Control under N<sub>2</sub> (g) at ~60 °C.
- h. The solution from (f.) or (g.) may be analyzed by color tests [see step (c.)]. A Negative Control will be analyzed along with the item extract(s).
- i. The solution from (f.) or (g.) may be analyzed by DART/TOFMS in positive and/or negative ionization mode (as appropriate based on the target analyte) by sampling the solution with the closed end of a glass capillary. Analyze the Negative Control(s), the Positive Control(s) (if applicable at this point), and PEG within the same data collection file. If the Positive Control(s) is determined after the initial DART/TOFMS analysis, then analyze the Positive Control(s) and PEG within a separate data collection file. It is also acceptable to analyze an item prior to extraction/dilution. For powder samples, a glass capillary can be wetted with deionized water and then touched to the sample; collect a blank glass capillary wetted with deionized water as a Negative Control.
- j. The solution from (f.) or (g.) may be analyzed by GC/MS in the electron impact (EI) mode. Also analyze the Negative Control and Positive Control(s) (if applicable at this point), and incorporate a solvent blank between each sample. If the Positive Control(s) is determined after the initial GC/MS analysis, then analyze the Positive Control(s) within a separate sequence (or edit the current sequence).
- k. The solution from (f.) or (g.) may be analyzed by GC/MS in the positive ion chemical ionization (PICI) or negative ion chemical ionization (NICI) mode as appropriate. Also analyze the Negative Control and a Positive Control (if applicable at this point), and incorporate a solvent blank between each sample. If the Positive Control(s) is determined after the initial GC/MS analysis, then analyze the Positive Control(s) within a separate sequence (or edit the current sequence). If an amphetamine (or other drug which gives a limited EI spectrum) is suspected, and FTIR or DART/TOFMS has not been utilized (or did not provide sufficient information), then PICI will be performed. A basic extraction using sodium bicarbonate and chloroform may improve the chromatography of methamphetamine (see Appendix A).

1. Other analytical techniques not listed above may be used to analyze a drug as deemed necessary {e.g., X-ray Powder Diffractometry (XRD), Liquid Chromatography/Mass Spectrometry [LC/MS, with electrospray (ESI) or atmospheric-pressure chemical ionization (APCI)], Raman spectrophotometry, Ultraviolet-Visible (UV-Vis) spectroscopy, Thin-Layer Chromatography (TLC)} provided that the instrumental conditions are retained in the case notes, the Negative Control and Positive Control samples provide the appropriate responses, and any relied upon references are retained in the case notes. These techniques are reserved for instances when the preceding steps don't provide sufficient data to identify a drug. If it is anticipated that the technique will be used routinely in the future for the drug, then the technique will be validated per the *Chemistry Unit Validation of Analytical Procedures* (CU QAOM 11).

## 7 Calculations

Not applicable.

## 8 Measurement Uncertainty

When quantitative results (e.g., weight, volume) are included in a *Laboratory Report*, measurement uncertainty will be estimated and reported following the *Chemistry Unit Procedures for Estimating Measurement Uncertainty* (CU QAOM 13). Uncertainty budget worksheets for each analytical balance approved for significant measurements are maintained electronically in CU.

## 9 Limitations

- The available sample size may limit or preclude some analytical techniques from being performed.
- Isomeric forms of a compound may not be differentiated by the techniques in this SOP. If relevant isomeric forms of a compound are not differentiated, this will be clearly stated in the *Laboratory Report*.

The following conclusions apply to the analysis of drugs:

- Identification (i.e. identified)
- Consistent with
- Not identified
- Inconclusive

Refer to *Chemistry Unit (CU) FBI Approved Standards for Scientific Testimony and Report Language for General Chemistry* (GenChem 32, ASSTR), *General Approach to Report Writing in General Chemistry* (GenChem 27), and *Department of Justice Uniform Language for Testimony and Reports for General Forensic Chemistry and Seized Drug Examinations* (GenChem ULTR) for examples of reporting examination conclusions and the associated limitations and decision criteria.

Refer to *General Chemistry Instrument Parameters* (GenChem 34) for instrumental limitations and decision criteria.

Refer to *General Chemistry Guidelines for Comparison of Mass Spectra* (GenChem 33) for mass spectra comparison decision criteria.

## 10 Safety

Take standard precautions for the handling of all chemicals, reagents, and standards. Some of the chemicals may be carcinogenic. Refer to the *FBI Laboratory Safety Manual* for the proper handling and disposal of all chemicals. Personal protective equipment should be used when handling any chemical and when performing any type of analysis.

## 11 References

Moffat AC, Osselton MD, Widdop B, Watts J. *Clarke's Analysis of Drugs and Poisons*, 4th ed., Pharmaceutical Press: 2011

Drug Enforcement Administration, Office of Forensic Sciences, *Analysis of Drugs Manual*, Revision 4, September 2019

*The Merck Index Online*, Royal Society of Chemistry

Scientific Working Group for the Analysis of Seized Drugs (SWGDRUG), *SWGDRUG Recommendations*, 8<sup>th</sup> Edition, 2019

Scientific Working Group for the Analysis of Seized Drugs (SWGDRUG), *SWGDRUG Monographs*, [www.swgdrug.org](http://www.swgdrug.org)

European Network of Forensic Science Institutes (ENFSI), *Guidelines on Representative Drug Sampling*, 2009



European Network of Forensic Science Institutes (ENFSI), *Guidelines on Sampling of Illicit Drugs for Qualitative Analysis*, 2<sup>nd</sup> Edition, 2016

*Analysis of Tablets and Capsules*; FBI Laboratory Chemistry Unit – General Chemistry SOP (GenChem 7)

*Sampling Guidelines for Bulk Materials and Multi-Unit Populations*; FBI Laboratory Chemistry Unit – General Chemistry SOP (GenChem 21)

*General Chemistry Instrument Parameters*; FBI Laboratory Chemistry Unit – General Chemistry SOP (GenChem 34)

*Guidelines for Comparison of Mass Spectra*; FBI Laboratory Chemistry Unit – General Chemistry SOP (GenChem 33)

*Chemistry Unit Procedures for Estimating Measurement Uncertainty*; FBI Laboratory Chemistry Unit – Quality Assurance and Operations Manual (CU QAOM 13)

*Chemistry Unit Validation of Analytical Procedures*; FBI Laboratory Chemistry Unit – Quality Assurance and Operations Manual (CU QAOM 11)

*FBI Laboratory Safety Manual*

Rev. #	Issue Date	History
2	01/15/20	Removed “Subunit” throughout. Removed previous section 1 (Introduction), section 3 (Principle), and section 7 (Calibration), and renumbered sections accordingly. Edited new section 1 for clarity and to include personnel. Defined ‘Chemistry Unit’ as ‘CU’. Changed lettered listing in section 2 to bullets and revised the list. Edited new sections 3.1 and 3.2 to add detail. Minor edits to sections 4 through 4.3 for clarity and to remove unnecessary detail. Section 5 edited to add more detail and to incorporate ‘non-statistical’ sampling language. Revised entirety of section 6 to include more detail and added DART/TOFMS as a routine technique. Minor edits made to section 8 for clarity. Changed the format of section 9 and added more detail to include flexibility for instrumental conditions. Revised section 10 for format and to include criteria for DART/TOFMS. Reformatted section 13 and updated content.
3	04/01/21	Section 1- added “and authorized” and edited to include only General Chemistry evidence. Section 2- removed ammonium hydroxide, added sodium sulfate (anhydrous). Section 3.2- minor grammatical edit (“into” to “within”). Section 5- added “on a heterogeneous item”. Section 6- added first sentence; step (a)- edited for clarity; step (c)- allowed for test tubes and other containers, as well as changes to reagent and sample amounts; step (f)- added last sentence (and added Appendix A); steps (j) and (k)- added “(or edit the current sequence)”; step (k)- added reference to Appendix A. Deleted previous sections 9 and 10. Section 8- edited last sentence to remove “s:\ drive”. Section 9- added the content below the first bulleted list. Section 11- added parenthetical details to CU references, added GenChem 34.

**Approval**

Redacted - Signatures on File

Chemistry Unit Chief:

Date: 03/31/2021

General Chemistry  
Technical Leader:

Date: 03/31/2021

## **Appendix A: Acid/Neutral and Alkaline Drug Extractions**

### **Acid/Neutral Drug Extraction-**

Mix several milligrams of a sample (homogenized if necessary) with several milliliters of deionized water in a test tube. Acidify the solution with 0.1 N hydrochloric acid until a pH of ~2 is achieved (check with pH paper). Add several milliliters of chloroform and rotate the mixture for approximately 10 minutes. Isolate the bottom chloroform layer and filter through pre-rinsed anhydrous sodium sulfate. Collect the chloroform layer into a labeled test tube and concentrate the solution under N<sub>2</sub> (g) flow at ~60 °C.

### **Alkaline Drug Extraction-**

Mix several milligrams of a sample (homogenized if necessary) with several milliliters of deionized water in a test tube. Add sodium bicarbonate until a pH of ~10 is achieved (check with pH paper). Add several milliliters of chloroform and rotate the mixture for approximately 10 minutes. Isolate the bottom chloroform layer and filter through pre-rinsed anhydrous sodium sulfate. Collect the chloroform layer into a labeled test tube and concentrate the solution under N<sub>2</sub> (g) flow at ~60 °C.