

Common Volatiles Analysis by Headspace GC-MSD/FID

Table of Contents

1	INTRODUCTION	3
2	SCOPE	3
3	PRINCIPLE	3
4	SPECIMEN CRITERIA	3
5	EQUIPMENT	3
5.1	Equipment	3
5.2	Consumables	4
5.3	Instruments	4
5.3.1	Columns	4
5.4	Software	4
5.5	Chemicals/Reagents	4
5.5.1	Purchased	4
5.5.2	Prepared	4
5.6	Standards/Controls	5
5.6.1	Purchased	5
5.6.2	Prepared	6
6	PROCEDURE	7
6.1	Screening/Identification by HS-GC/MSD	7
6.2	Confirmation/Quantitation by HS-GC/FID	8
6.3	Screening or Confirmation for Nonstandard Samples	9
7	ANALYTICAL PARAMETERS	9
7.1	Pipettor Settings (MSD and FID)	9
7.2	Mass Spectrometry Method (Screening)	9
7.2.1	Inlets	9
7.2.2	Columns	10
7.2.3	Oven	11
7.2.4	Detector	11
7.2.5	Aux Heaters	11
7.2.6	Column Configuration	11
7.2.7	Module Configuration	12
7.2.8	GC Readiness	12
7.2.9	Quadrupole Settings	12
7.3	Gerstel AutoSampler Settings (FID and MSD)	13
7.4	Flame Ionization Method (Confirmation)	14
7.4.1	Inlets	14
7.4.2	Columns	14

7.4.3	Oven	15
7.4.4	Detector	15
7.4.5	Aux Heater	15
7.4.6	Signals	16
7.4.7	Column Configuration.....	16
7.4.8	Module Configuration.....	16
7.4.9	GC Readiness.....	17
8	DATA ANALYSIS	18
8.1	Decision Criteria.....	18
8.1.1	FID Method	18
8.1.2	MSD Method.....	19
8.1.3	Batch Acceptance.....	20
8.2	Calculations.....	20
8.2.1	MSD Screening	20
8.2.2	FID Confirmation.....	20
8.2.3	Characterization of Whole Blood Controls.....	20
9	REPORTING.....	20
9.1	Measurement Uncertainty	20
9.2	MSD Screening	21
9.3	FID Confirmation.....	21
9.4	Quantitative Values	21
10	CORRECTIVE MEASURES.....	21
11	PERFORMANCE CHARACTERISTICS	21
11.1	LOD (MSD and FID)	21
11.2	LOQ (FID).....	22
11.3	Linearity (FID).....	22
11.4	Precision (FID)	22
11.5	Carryover.....	22
11.6	Processed Sample Stability	22
12	LIMITATIONS	22
12.1	Interferences.....	22
12.2	Interpretation	22
13	SAFETY	23
14	REVISION HISTORY	23

Common Volatiles Analysis by Headspace GC-MSD/FID

1 INTRODUCTION

This procedure allows for the screening, identification, confirmation, and quantitation of common volatile chemicals.

2 SCOPE

Analyses	<input checked="" type="checkbox"/> Screening <input checked="" type="checkbox"/> Confirmation <input checked="" type="checkbox"/> Quantitation
Matrices	Blood, serum, urine, vitreous fluid and other liquids
Analytes	Ethanol, acetone, isopropanol, methanol (Target Compounds)
Personnel	This document applies to authorized personnel who perform the described tasks, singly or in combination.

3 PRINCIPLE

Sample and a diluent containing internal standard are added to a headspace vial using a pipette operating in dilute mode. Samples are qualitatively screened for target compounds by headspace gas chromatography with mass spectrometry (HS-GC/MS). Specimens are quantitatively confirmed through a separate analysis using headspace gas chromatography with flame ionization detection (HS-GC/FID). The headspace technique is based on Henry's gas law, which states that, for dilute solutions at a constant temperature and at equilibrium, a direct relationship exists between the amount of volatile analyte dissolved in a liquid and the amount of the analyte in the vapor above the solution. Determination of the concentration of the volatile analyte in the sampled headspace allows for the calculation of the concentration of that analyte in the original liquid sample.

4 SPECIMEN CRITERIA

0.1 mL of sample is used per replicate analysis.

5 EQUIPMENT

Use of equivalent equipment is allowable.

5.1 Equipment

Electronic Pipettor	Eppendorf Xplorer Plus, single channel, 50-1000 μ L range
Headspace vial cap crimper	Standard, 10mL
Routine Laboratory Glassware and supplies	Volumetric flasks (50, 100 and 1000 mL), pipettes, disposable tissue grinder
Laboratory Balance	Standard, \geq 0.1g resolution. Traceable.

5.2 Consumables

Inlet Liner	Restek 1.0mm Topaz Straight Liner. PN 23333
Pipette Tips	Biotix™ uTIP™ Filter Pipette Tips for Universal Pipettes, Standard. PN M12509FC96
Headspace vials	Gerstel, crimp cap vials, 10 mL, 100 pack. PN 093640-005-00.
Headspace vial caps, magnetic	Gerstel, crimp caps with septum for vials, 100 pack. PN 093640-063-00
Storage tubes	Thermo Matrix 1.0mL ScrewTop tubes PN 3741WHI

5.3 Instruments

GC/MS with Headspace Autosampler	El ionization, Gerstel autosampler
----------------------------------	------------------------------------

5.3.1 Columns

GC Columns	Restek RTX-BAC Plus 1: 30m X 0.32mm X 1.8 μm PN 18004 Restek RTX-BAC Plus 2: 30m X 0.32mm X 0.6 μm PN 18006
------------	--

5.4 Software

Component	Software	Version
Operating System	Microsoft Windows	10 Pro
GC/MS	Agilent 7890B	B.02.05.2
	MassHunter Workstation	10.0.368
Autosampler	Gerstel Maestro 1	1.5.3.74/3.5
Data Analysis	MassHunter Workstation Quantitative Analysis	10.0

5.5 Chemicals/Reagents

5.5.1 Purchased

Item	Supplier*	Description**	Part Number*
t-butanol	Sigma-Aldrich	ACS/Reagent Grade	360538
Ethanol	Sigma-Aldrich	HPLC grade	E7023
Methanol	Fisher Scientific	HPLC grade	A454
Isopropanol	Fisher Scientific	HPLC grade	A451
Acetone	Fisher Scientific	HPLC grade	A949
Deionized water	Laboratory supplied	18.2 MΩ•cm	N/A

*use of an equivalent product is allowable **listed grade or better

5.5.2 Prepared

Depending upon the batch size, the absolute amounts may be adjusted so long as the ratios of components are maintained.

5.5.2.1 Stock Sample Diluent (1.0 g/100mL)

Step	Action	Amount	Component/Information
1	Acquire	1	Volumetric flask, 100 mL
2	Add	~90 mL	Deionized water
3	Add	1.0 g	t-butanol
4	QS	100 mL	Deionized water
5	Mix		
6	Transfer		Glass container
7	Storage		Refrigerated or ambient
8	Stability		≥ 6 months
9	Prepares	100 mL	(20 Sample Diluent preparations)

5.5.2.2 Sample Diluent (0.005 g/100mL)

Step	Action	Amount	Component/Information
1	Acquire	1	Volumetric flask, 1000 mL
2	Add	5.0 mL	Stock Sample Diluent
3	QS	1000 mL	Deionized water
4	Mix		
5	Transfer		Glass container, tightly sealed
6	Storage		Ambient. Smaller satellite container may also be used.
7	Stability		≥ 6 months
8	Prepares	1000 mL	(1250 analyses)

5.6 Standards/Controls

5.6.1 Purchased

Analyte	Supplier *	Description	Part Number*
Multicomponent Volatiles	Cerilliant	C1-C6 levels containing ethanol, methanol, isopropanol and acetone at 0.010, 0.025, 0.050, 0.100, 0.200 and 0.400 g/100mL in water	A-127
Ethanol	Cerilliant	C7 level containing ethanol at 0.500 g/100mL in water	E-053
Multicomponent Volatiles	Cliniqa	Contain ethanol, methanol, isopropanol, and acetone in whole human blood (varying concentrations)	93221, 93222

*Use of an equivalent product is allowable; be sure to consider additional validation/verification and update Measurement Traceability records if alternate calibrators are used. Store refrigerated. Stability determined by manufacturer. Cerilliant solutions will be used shortly after opening and not reused on subsequent days. Cliniqa controls will be used within 45 days of opening. Refer to 5.5 for additional starting materials.

5.6.2 Prepared

Depending upon the batch size, the absolute amounts may be adjusted so long as the ratios of components are maintained.

At first opening of each new Cliniqua control vial, the contents of the vial will be portioned out into storage tubes that limit the headspace during storage. The date of opening and portioning will be recorded on each storage vial. Each of these vials should be used on the day of opening and then discarded.

5.6.2.1 TOX200 Stock System Suitability Sample (0.100 g/100mL)

Step	Action	Amount	Component/Information
1	Acquire	1	Volumetric flask, 50 mL
2	Add	~25 mL	Deionized Water
3	Add	0.064 mL	Each of stock ethanol, acetone, isopropanol, methanol
3	QS	50 mL	Deionized water
4	Mix		
5	Transfer		Glass container, tightly sealed
6	Storage		Refrigerated.
7	Stability		≥ 12 months
8	Prepares	50 mL	Of stock material

5.6.2.2 TOX200 System Suitability Sample (0.010 g/100mL, S³)



Step	Action	Amount	Component/Information
1	Acquire	1	Volumetric flask, 50 mL
2	Add	~25 mL	Deionized Water
3	Add	5 mL	Stock System Suitability Sample
3	QS	50 mL	Deionized water
4	Mix		
5	Transfer		Glass container, tightly sealed
6	Storage		Refrigerated.
7	Stability		≥ 12 months
8	Prepares	50 mL	(500 analyses)

Analysis of an S³ is used to verify system performance for both FID and MSD methods prior to case analysis.

6 PROCEDURE



6.1 Screening/Identification by HS-GC/MSD

Batch Building: Samples used for screening may be aliquoted in advance up to 15 days prior to analysis (samples and diluent are portioned into a headspace vial, sealed, and placed in secure refrigerated storage). The same lot of Sample Diluent is used throughout a given batch. Any calibrated Xplorer Plus pipette may be used. Batches are coded according to the scheme TOX200.YYYYMMDD-MSD.

Step		Activity	Note	Reference/Lot																																																																																
6.1.1	<input type="checkbox"/>	<p>Samples: Using an Eppendorf Xplorer pipette fitted with a tip, aliquot 800 µL of Sample Diluent and 100 µL of sample into a 10 mL headspace vial. Crimp vial firmly using a magnetic cap. Use a new tip for each sampling.</p>	See pipette settings in Section 7.1	<p>S³ Sample Diluent Xplorer Plus Pipette</p> 																																																																																
6.1.2	<input type="checkbox"/>	<p>Quality Control Materials: To start a batch, pipet the following QC materials:</p> <ul style="list-style-type: none"> Negative Control (deionized water) 0.010 g% (CRM) 0.200 g% (CRM) <p>Upon aliquot of the final case sample for the batch, include a closing control:</p> <ul style="list-style-type: none"> 0.100 g% (CRM) 		<p>Negative Control Calibrator C1-C6 Set</p> 																																																																																
6.1.3	<input type="checkbox"/>	<p>Batch Analysis: Input the samples into the instrument sequence using the following order and format:</p> <table border="1"> <thead> <tr> <th>Vial</th> <th>Sample Type</th> <th>Sample Name</th> <th>Method File</th> <th>Data Path</th> <th>Data File</th> <th>Tray</th> <th>Volume</th> </tr> </thead> <tbody> <tr> <td>1</td> <td>Sample</td> <td>NEG</td> <td>TOX200-MSD.M</td> <td>D:\MassHunter\Data\TOX200\20201019</td> <td>20201019MSD-01</td> <td>Rack 1,R60/10-CVM</td> <td>1000.0</td> </tr> <tr> <td>2</td> <td>Cal</td> <td>CAL 0.010</td> <td>TOX200-MSD.M</td> <td>D:\MassHunter\Data\TOX200\20201019</td> <td>20201019MSD-02</td> <td>Rack 1,R60/10-CVM</td> <td>1000.0</td> </tr> <tr> <td>3</td> <td>Cal</td> <td>CAL 0.200</td> <td>TOX200-MSD.M</td> <td>D:\MassHunter\Data\TOX200\20201019</td> <td>20201019MSD-03</td> <td>Rack 1,R60/10-CVM</td> <td>1000.0</td> </tr> <tr> <td>4</td> <td>Sample</td> <td>Case Sample 1</td> <td>TOX200-MSD.M</td> <td>D:\MassHunter\Data\TOX200\20201019</td> <td>20201019MSD-04</td> <td>Rack 1,R60/10-CVM</td> <td>1000.0</td> </tr> <tr> <td>5</td> <td>Sample</td> <td>Case Sample 2</td> <td>TOX200-MSD.M</td> <td>D:\MassHunter\Data\TOX200\20201019</td> <td>20201019MSD-05</td> <td>Rack 1,R60/10-CVM</td> <td>1000.0</td> </tr> <tr> <td>6</td> <td>Sample</td> <td>Case Sample 3</td> <td>TOX200-MSD.M</td> <td>D:\MassHunter\Data\TOX200\20201019</td> <td>20201019MSD-06</td> <td>Rack 1,R60/10-CVM</td> <td>1000.0</td> </tr> <tr> <td>7</td> <td>Sample</td> <td>Case Sample 4</td> <td>TOX200-MSD.M</td> <td>D:\MassHunter\Data\TOX200\20201019</td> <td>20201019MSD-07</td> <td>Rack 1,R60/10-CVM</td> <td>1000.0</td> </tr> <tr> <td>8</td> <td>Sample</td> <td>Case Sample 5...</td> <td>TOX200-MSD.M</td> <td>D:\MassHunter\Data\TOX200\20201019</td> <td>20201019MSD-08</td> <td>Rack 1,R60/10-CVM</td> <td>1000.0</td> </tr> <tr> <td>9</td> <td>QC</td> <td>CONTROL 0.10</td> <td>TOX200-MSD.M</td> <td>D:\MassHunter\Data\TOX200\20201019</td> <td>20201019MSD-09</td> <td>Rack 1,R60/10-CVM</td> <td>1000.0</td> </tr> </tbody> </table> <p>A maximum of 116 samples may be analyzed in one batch. In batches that contain 20 or more samples, an additional “mid-run” positive control will be added.</p>	Vial	Sample Type	Sample Name	Method File	Data Path	Data File	Tray	Volume	1	Sample	NEG	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-01	Rack 1,R60/10-CVM	1000.0	2	Cal	CAL 0.010	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-02	Rack 1,R60/10-CVM	1000.0	3	Cal	CAL 0.200	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-03	Rack 1,R60/10-CVM	1000.0	4	Sample	Case Sample 1	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-04	Rack 1,R60/10-CVM	1000.0	5	Sample	Case Sample 2	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-05	Rack 1,R60/10-CVM	1000.0	6	Sample	Case Sample 3	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-06	Rack 1,R60/10-CVM	1000.0	7	Sample	Case Sample 4	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-07	Rack 1,R60/10-CVM	1000.0	8	Sample	Case Sample 5...	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-08	Rack 1,R60/10-CVM	1000.0	9	QC	CONTROL 0.10	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-09	Rack 1,R60/10-CVM	1000.0		
Vial	Sample Type	Sample Name	Method File	Data Path	Data File	Tray	Volume																																																																													
1	Sample	NEG	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-01	Rack 1,R60/10-CVM	1000.0																																																																													
2	Cal	CAL 0.010	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-02	Rack 1,R60/10-CVM	1000.0																																																																													
3	Cal	CAL 0.200	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-03	Rack 1,R60/10-CVM	1000.0																																																																													
4	Sample	Case Sample 1	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-04	Rack 1,R60/10-CVM	1000.0																																																																													
5	Sample	Case Sample 2	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-05	Rack 1,R60/10-CVM	1000.0																																																																													
6	Sample	Case Sample 3	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-06	Rack 1,R60/10-CVM	1000.0																																																																													
7	Sample	Case Sample 4	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-07	Rack 1,R60/10-CVM	1000.0																																																																													
8	Sample	Case Sample 5...	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-08	Rack 1,R60/10-CVM	1000.0																																																																													
9	QC	CONTROL 0.10	TOX200-MSD.M	D:\MassHunter\Data\TOX200\20201019	20201019MSD-09	Rack 1,R60/10-CVM	1000.0																																																																													

6.2 Confirmation/Quantitation by HS-GC/FID

Batch Building: Allow specimens and quality control samples to stand at room temperature for at least 15 minutes. Samples used for confirmation/quantitation are aliquoted from the original item or an intermediary container into a headspace vial and sealed. The same lot of Sample Diluent is used throughout a given batch. Any calibrated Xplorer Plus pipette may be used. Batches are coded according to the scheme TOX200.YYYYMMDD-FID.


Step		Activity	Note	Reference/Lot																																																																																																																																																																								
6.2.1	<input type="checkbox"/>	<p>Samples: Using an Eppendorf Xplorer pipette fitted with a filter tip, aliquot 800 µL of Sample Diluent and 100 µL of sample into a 10 mL headspace vial. Crimp vial firmly using a magnetic cap. Perform in duplicate. Use a new tip for each sampling.</p>	See pipette settings in Section 7.1	<p>S³ Sample Diluent Xplorer Plus Pipette</p> 																																																																																																																																																																								
6.2.2	<input type="checkbox"/>	<p>Quality Control Materials: Use the following QC materials for each batch:</p> <ul style="list-style-type: none"> Negative Control (deionized water) CAL1-CAL6, CAL7 (CRM) Cliniqa Controls (Two Levels) 		<p>Negative Control Calibrator C1-C6 Set Calibrator C7 Cliniqa Level 1 Cliniqa Level 2</p> 																																																																																																																																																																								
6.2.3	<input type="checkbox"/>	<p>Batch Analysis: Input the samples into the instrument sequence using the following order and format:</p> <table border="1"> <thead> <tr> <th>Vial</th> <th>Sample Type</th> <th>Sample Name</th> <th>Method File</th> <th>Data Path</th> <th>Data File</th> <th>Tray</th> <th>Volume</th> </tr> </thead> <tbody> <tr><td>1</td><td>Sample</td><td>Negative QC</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-01</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>2</td><td>Cal</td><td>CAL1</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-02</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>3</td><td>Cal</td><td>CAL2</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-03</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>4</td><td>Cal</td><td>CAL3</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-04</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>5</td><td>Cal</td><td>CAL4</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-05</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>6</td><td>Cal</td><td>CAL5</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-06</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>7</td><td>Cal</td><td>CAL6</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-07</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>8</td><td>Cal</td><td>CAL7</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-08</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>9</td><td>Sample</td><td>blank</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-09</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>10</td><td>QC</td><td>Low QC</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-10</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>11</td><td>QC</td><td>High QC</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-11</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>12</td><td>Sample</td><td>blank</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-12</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>13</td><td>Sample</td><td>Case 1</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-13</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>14</td><td>Sample</td><td>Case 1</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-14</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>15</td><td>Sample</td><td>blank</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-15</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>16</td><td>Sample</td><td>Case 2</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-16</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>17</td><td>Sample</td><td>Case 2...</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-17</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>18</td><td>Sample</td><td>blank</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-18</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>19</td><td>QC</td><td>Low QC</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-19</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> <tr><td>20</td><td>QC</td><td>High QC</td><td>TOX200-FID.M</td><td>D:\MassHunter\Data\TOX200\20201019</td><td>20201019FID-20</td><td>Rack 1,R60/10-CVM</td><td>1000.0</td></tr> </tbody> </table> <p>A maximum of 35 samples may be analyzed in one batch. In batches that contain 20 or more samples, an additional “mid-run” positive control will be added.</p>	Vial	Sample Type	Sample Name	Method File	Data Path	Data File	Tray	Volume	1	Sample	Negative QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-01	Rack 1,R60/10-CVM	1000.0	2	Cal	CAL1	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-02	Rack 1,R60/10-CVM	1000.0	3	Cal	CAL2	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-03	Rack 1,R60/10-CVM	1000.0	4	Cal	CAL3	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-04	Rack 1,R60/10-CVM	1000.0	5	Cal	CAL4	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-05	Rack 1,R60/10-CVM	1000.0	6	Cal	CAL5	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-06	Rack 1,R60/10-CVM	1000.0	7	Cal	CAL6	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-07	Rack 1,R60/10-CVM	1000.0	8	Cal	CAL7	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-08	Rack 1,R60/10-CVM	1000.0	9	Sample	blank	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-09	Rack 1,R60/10-CVM	1000.0	10	QC	Low QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-10	Rack 1,R60/10-CVM	1000.0	11	QC	High QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-11	Rack 1,R60/10-CVM	1000.0	12	Sample	blank	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-12	Rack 1,R60/10-CVM	1000.0	13	Sample	Case 1	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-13	Rack 1,R60/10-CVM	1000.0	14	Sample	Case 1	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-14	Rack 1,R60/10-CVM	1000.0	15	Sample	blank	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-15	Rack 1,R60/10-CVM	1000.0	16	Sample	Case 2	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-16	Rack 1,R60/10-CVM	1000.0	17	Sample	Case 2...	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-17	Rack 1,R60/10-CVM	1000.0	18	Sample	blank	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-18	Rack 1,R60/10-CVM	1000.0	19	QC	Low QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-19	Rack 1,R60/10-CVM	1000.0	20	QC	High QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-20	Rack 1,R60/10-CVM	1000.0		
Vial	Sample Type	Sample Name	Method File	Data Path	Data File	Tray	Volume																																																																																																																																																																					
1	Sample	Negative QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-01	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
2	Cal	CAL1	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-02	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
3	Cal	CAL2	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-03	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
4	Cal	CAL3	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-04	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
5	Cal	CAL4	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-05	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
6	Cal	CAL5	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-06	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
7	Cal	CAL6	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-07	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
8	Cal	CAL7	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-08	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
9	Sample	blank	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-09	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
10	QC	Low QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-10	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
11	QC	High QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-11	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
12	Sample	blank	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-12	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
13	Sample	Case 1	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-13	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
14	Sample	Case 1	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-14	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
15	Sample	blank	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-15	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
16	Sample	Case 2	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-16	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
17	Sample	Case 2...	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-17	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
18	Sample	blank	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-18	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
19	QC	Low QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-19	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					
20	QC	High QC	TOX200-FID.M	D:\MassHunter\Data\TOX200\20201019	20201019FID-20	Rack 1,R60/10-CVM	1000.0																																																																																																																																																																					

6.3 Screening or Confirmation for Nonstandard Samples

A sample is considered nonstandard if it cannot be rendered homogenous through mixing/vortexing, which are the preferred methods. If a case sample is clotted and cannot be pipetted accurately, it may be homogenized with a clot grinder before pipetting. If the entire specimen or a portion of the specimen is homogenized, the details of the homogenization are recorded. For additional guidance, refer to TOX-100.

7 ANALYTICAL PARAMETERS

7.1 Pipettor Settings (MSD and FID)



7.2 Mass Spectrometry Method (Screening)

7.2.1 Inlets

Select...	Split-Splitless Inlet	Select Liner...	Liner: Restek 23313; Lot 111209-1: 500 µL (Topaz Splitless straight liner)															
ALS																		
▼ Inlets																		
SSL - Front																		
SSL - Back																		
Columns																		
Oven																		
Detectors																		
Aux Heaters																		
Events																		
Signals																		
▼ Configuration																		
Miscellaneous																		
Columns																		
Modules																		
ALS																		
Backflush																		
Readiness																		
GC Calculators																		
	<table border="1"><thead><tr><th></th><th>Actual</th><th>Setpoint</th></tr></thead><tbody><tr><td><input checked="" type="checkbox"/> Heater:</td><td>200 °C</td><td>200 °C</td></tr><tr><td><input checked="" type="checkbox"/> Pressure:</td><td>1.434 psi</td><td>1.434 psi</td></tr><tr><td>Total Flow:</td><td>18.62 mL/min</td><td>18.623 mL/min</td></tr><tr><td><input checked="" type="checkbox"/> Septum Purge Flow:</td><td>3 mL/min</td><td>3 mL/min</td></tr></tbody></table>		Actual	Setpoint	<input checked="" type="checkbox"/> Heater:	200 °C	200 °C	<input checked="" type="checkbox"/> Pressure:	1.434 psi	1.434 psi	Total Flow:	18.62 mL/min	18.623 mL/min	<input checked="" type="checkbox"/> Septum Purge Flow:	3 mL/min	3 mL/min		
	Actual	Setpoint																
<input checked="" type="checkbox"/> Heater:	200 °C	200 °C																
<input checked="" type="checkbox"/> Pressure:	1.434 psi	1.434 psi																
Total Flow:	18.62 mL/min	18.623 mL/min																
<input checked="" type="checkbox"/> Septum Purge Flow:	3 mL/min	3 mL/min																
	<table border="1"><tr><td>▲ Inlet Mode (Split 10 : 1)</td><td></td><td></td></tr><tr><td>Split</td><td>Split Ratio:</td><td></td></tr><tr><td></td><td>10 : 1</td><td>Split Flow 14.203 mL/min</td></tr></table>	▲ Inlet Mode (Split 10 : 1)			Split	Split Ratio:			10 : 1	Split Flow 14.203 mL/min								
▲ Inlet Mode (Split 10 : 1)																		
Split	Split Ratio:																	
	10 : 1	Split Flow 14.203 mL/min																
	<table border="1"><tr><td>▲ Gas Saver (Off)</td><td></td><td></td></tr><tr><td><input type="checkbox"/> On</td><td>20 mL/min</td><td>After: 2 min</td></tr></table>	▲ Gas Saver (Off)			<input type="checkbox"/> On	20 mL/min	After: 2 min											
▲ Gas Saver (Off)																		
<input type="checkbox"/> On	20 mL/min	After: 2 min																

Select... **Split-Splitless Inlet** Select Liner... Liner: Restek 23313; Lot 111209-1: 500 µL (Topaz Splitless straight liner)

	Actual	Setpoint
<input checked="" type="checkbox"/> Heater:	200 °C	200 °C
<input checked="" type="checkbox"/> Pressure:	15.748 psi	15.749 psi
Total Flow:	47 mL/min	47 mL/min
<input checked="" type="checkbox"/> Septum Purge Flow:	3.001 mL/min	3 mL/min

Inlet Mode (Split 10 : 1)

Split Ratio: 10 : 1 Split Flow 40 mL/min

Gas Saver (Off)

On 20 mL/min After: 2 min

7.2.2 Columns

Select...	#	Selection	Columns
ALS		Front SS Inlet He ---> Restek 18004: 2037	Columns Control Mode <input checked="" type="checkbox"/> On Flow: 1.42 mL/min (Actual) / 1.4203 mL/min (Setpoint) Pressure: 1.433 psi / 1.434 psi Average Velocity: 43.215 cm/sec Holdup Time: 1.157 min Constant Flow Post Run: 0.44594 mL/min Column #1 Configuration Change Column... Calibrate Column... Lock Column...
Inlets SSL - Front SSL - Back Columns	1	RTX-BAC Plus 1 5 °C—240 °C (260 °C): 30 m x 320 µm x 1.8 µm ---> MSD	
Oven	2	Back SS Inlet He ---> Restek 18006: 2038	((Initial): 0 min He @ 40 °C Oven Out: MSD 30 m x 320 µm x 1.8 µm
Detectors Aux Heaters Events Signals Configuration Miscellaneous Columns Modules ALS Backflush Readiness		RTX-BAC Plus 2 5 °C—240 °C (260 °C): 30 m x 320 µm x 0.6 µm ---> Front Detector FID	
		Aux EPC 1 He	
		Aux EPC 2 He	
		Aux EPC 3 He	

Select...	#	Selection	Columns
ALS		Front SS Inlet He ---> Restek 18004: 2037	Columns Control Mode <input checked="" type="checkbox"/> On Flow: 4 mL/min (Actual) / 4 mL/min (Setpoint) Pressure: 15.748 psi / 15.749 psi Average Velocity: 54.889 cm/sec Holdup Time: 0.91093 min Constant Flow Post Run: 0.71303 mL/min Column #2 Configuration Change Column... Calibrate Column... Lock Column...
Inlets SSL - Front SSL - Back Columns	1	RTX-BAC Plus 1 5 °C—240 °C (260 °C): 30 m x 320 µm x 1.8 µm ---> MSD	
Oven	2	Back SS Inlet He ---> Restek 18006: 2038	((Initial): 0 min He @ 40 °C Oven Out: Ambient Pressure 30 m x 320 µm x 0.6 µm
Detectors Aux Heaters Events Signals Configuration Miscellaneous Columns Modules ALS Backflush Readiness		RTX-BAC Plus 2 5 °C—240 °C (260 °C): 30 m x 320 µm x 0.6 µm ---> Front Detector FID	
		Aux EPC 1 He	
		Aux EPC 2 He	
		Aux EPC 3 He	

7.2.3 Oven

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven**
- Detectors
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous
 - Columns

Actual

Oven Temp On

40 °C

Equilibration Time

Maximum Oven Temperature

Override Column Max: 260 °C

	Rate °C/min	Value °C	Hold Time min	Run Time min
▶ (Initial)		40	4	4
*				

Post Run:

Post Run Time:

7.2.4 Detector

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors**
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous

FID

	Actual	Setpoint
<input checked="" type="checkbox"/> Heater:	250 °C	<input style="width: 60px;" type="text" value="250 °C"/>
<input type="checkbox"/> Air Flow:	-0.4557 mL/min	<input style="width: 60px;" type="text" value="400 mL/min"/>
<input type="checkbox"/> H2 Fuel Flow:	0.1125 mL/min	<input style="width: 60px;" type="text" value="30 mL/min"/>
<input type="checkbox"/> Makeup Flow: (He)	-0.05167 mL/min	<input style="width: 60px;" type="text" value="25 mL/min"/>
▶ Carrier Gas Flow Correction (None)		
<input type="checkbox"/> Flame	0 pA	
▶ No Column Comp		

7.2.5 Aux Heaters

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors
- Aux Heaters**

Aux Heaters

Thermal Aux 2 (MSD Transfer Line)

Actual

On

250 °C

7.2.6 Column Configuration

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous
 - Columns**

Flow Paths:
Front Inlet---->#1---->MSD
Back Inlet---->#2---->Front Detector

Options

Column Outlet Pressure:

	Column	Calibration Results	Inlet	Outlet	Heated By
↑	1 Restek 18004: 2037 RTX-BAC Plus 1 5 °C—240 °C (260 °C): 30 m x 320 µm x 1.8 µm	Uncalibrated	Front Inlet	MSD	Oven
↓	2 Restek 18006: 2038 RTX-BAC Plus 2 5 °C—240 °C (260 °C): 30 m x 320 µm x 0.6 µm	Uncalibrated	Back Inlet	Front Detector	Oven
	3 No Column Installed	Uncalibrated	Unspecified	Other	Oven
	4 No Column Installed	Uncalibrated	Unspecified	Other	Oven

7.2.7 Module Configuration

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous
 - Columns
 - Modules**
 - ALS
 - Backflush
 - Readiness
 - GC Calculators

Front Inlet
SS Inlet He

Back Inlet
SS Inlet He

Front Detector
FID
Makeup N2
Set Lit Offset with GC Keyboard.

Aux EPC 1,2,3
Aux EPC 1 He

Aux EPC 1,2,3
Aux EPC 2 He

Aux EPC 1,2,3
Aux EPC 3 He

7.2.8 GC Readiness

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors
- Aux Heaters
- Events
- Signals

GC Readiness

Only checked components will affect the GC readiness

- Oven
- Front Inlet (SS Inlet)
- Back Inlet (SS Inlet)
- Front Detector (FID)
- Aux EPC 1
- Aux EPC 2
- Aux EPC 3
- Thermal Aux 2 (MSD Transfer Line)

7.2.9 Quadrupole Settings

Single Quadrupole MS Method Editor

Tune File: lomass.u

Tune Type: EI
Tune EMV: 901
CI Gas Valve: ---
CI Flow: --- %

MS Source: 230 (Actual) 230 (Setpoint)
MS Quad: 150 (Actual) 150 (Setpoint) **Apply**

Acquisition Type: Scan

Run Time: 10.00 min
Solvent Delay: 1.40 min

Detector Setting
 Trace Ion Detection

EM Setting: Gain Factor
Gain Factor: 1.000
Applied EM Voltage (V): 930
 EM Saver
Limit: Sum Limit 1e8 (Default)

Scan Time Segments

	Time	Start Mass	End Mass	Threshold	Scan Speed (u/s)	Frequency (scans/sec)	Cycle Time (ms)	Step Size (m/z)
▶	1.40	27.00	100.00	150	1,562 [N=2]	14.4	69.40	0.1

SIM Time Segments

	Time	Group Name	Number of Ions	Total Dwell Time (ms)	Cycle Time (Hz)	Resolution	Gain Factor	Calculated EMV
▶	1.40		1	1	100	8.3333	Low	

SIM Real-Time Plot Timed Events

	m/z	Dwell Time	Plot Ion	Label
▶	74.10	100	<input checked="" type="checkbox"/>	
*			<input type="checkbox"/>	

7.3 Gerstel AutoSampler Settings (FID and MSD)

Settings

Runtime min

Minimum Runtime: 4.00 min (limited by GC)

GC Cool Down Time min

Time required to cool down the GC Oven to initial temperature

Settings

Runtime min

Minimum Runtime: 3.00 min (limited by GC)

GC Cool Down Time min

Time required to cool down the GC Oven to initial temperature

MSD (Screening)

FID (Confirmation)

Headspace Injection Settings Options

Syringe Settings

Syringe ▾

Syringe Temp. (°C)

Flush Time (s)

Sample

Inj. Volume (µL)

Inj. Speed (µL/s)

Pullup Delay (s)

Fill Volume (µL)

Fill Strokes

Fill Speed (µL/s)

Pre Inj. Delay (s)

Post Inj. Delay (s)

Inj. Penetration (mm)

Sample Tray Type ▾

Vial Penetration (mm)

Sample Preparation

Sample Mode ▾

— Heating and Incubation —

Incubator ▾

Incubation Temp. (°C) **50**

Incubation Time (min)

Agitator On Time (s)

Agitator Off Time (s)

Agitator Speed (rpm)

Headspace Injection Settings Options

Multiple Headspace Sample Enrichment (MHSE) and/or Pressurize

Pressurize Sample

Injections per Run

Delay Time (min)

7.4 Flame Ionization Method (Confirmation)

7.4.1 Inlets

Select... **Split-Splitless Inlet** Liner: Restek 23313; Lot 111209-1: 500 µL (Topaz Splitless straight liner)

ALS
 Inlets
 SSL - Front
 SSL - Back
 Columns
 Oven
 Detectors
 Aux Heaters
 Events
 Signals
 Configuration
 Miscellaneous
 Columns
 Modules
 ALS
 Backflush
 Readiness
 GC Calculators

	Actual	Setpoint
<input checked="" type="checkbox"/> Heater:	200 °C	200 °C
<input checked="" type="checkbox"/> Pressure:	1.433 psi	1.4338 psi
Total Flow:	6.42 mL/min	6.4203 mL/min
<input checked="" type="checkbox"/> Septum Purge Flow:	2 mL/min	2 mL/min

Septum Purge Flow Mode: Standard

Inlet Mode (Splitless)
 Splitless

Gas Saver (Off)
 On 20 mL/min After: 3 min

Select... **Split-Splitless Inlet** Liner: Restek 23313; Lot 111209-1: 500 µL (Topaz Splitless straight liner)

ALS
 Inlets
 SSL - Front
 SSL - Back
 Columns
 Oven
 Detectors
 Aux Heaters
 Events
 Signals
 Configuration
 Miscellaneous
 Columns
 Modules
 ALS
 Backflush
 Readiness
 GC Calculators

	Actual	Setpoint
<input checked="" type="checkbox"/> Heater:	200 °C	200 °C
<input checked="" type="checkbox"/> Pressure:	15.748 psi	15.749 psi
Total Flow:	47 mL/min	47 mL/min
<input checked="" type="checkbox"/> Septum Purge Flow:	3 mL/min	3 mL/min

Inlet Mode (Split 10 : 1)
 Split

Gas Saver (Off)
 On 20 mL/min After: 1 min

7.4.2 Columns

Select... **Columns**

#	Selection
1	Front SS Inlet He ---> Restek 18004: 2037 RTX-BAC Plus 1 5 °C—240 °C (260 °C): 30 m x 320 µm x 1.8 µm ---> MSD
2	Back SS Inlet He ---> Restek 18006: 2038 RTX-BAC Plus 2 5 °C—240 °C (260 °C): 30 m x 320 µm x 0.6 µm ---> Front Detector FID
	Aux EPC 1 He
	Aux EPC 2 He
	Aux EPC 3 He

ALS
 Inlets
 SSL - Front
 SSL - Back
 Columns
 Oven
 Detectors
 Aux Heaters
 Events
 Signals
 Configuration
 Miscellaneous
 Columns
 Modules
 ALS
 Backflush
 Readiness
 GC Calculators

Columns

Control Mode
 On

	Actual	Setpoint
Flow	1.42 mL/min	1.4203 mL/min
Pressure	1.434 psi	1.4338 psi
Average Velocity		43.215 cm/sec
Holdup Time		1.157 min

(Initial): 0 min
 He @ 40 °C Oven
 Out: MSD
 30 m x 320 µm x 1.8 µm

Constant Flow

Post Run: 0.44594 mL/min

Column #1 Configuration

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns**
- Oven
- Detectors
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous
 - Columns
 - Modules
 - ALS
 - Backflush
 - Readiness
 - GC Calculators

#	Selection
1	Front SS Inlet He ---> Restek 18004: 2037 RTX-BAC Plus 1 5 °C—240 °C (260 °C); 30 m x 320 µm x 1.8 µm ---> MSD
2	Back SS Inlet He ---> Restek 18006: 2038 RTX-BAC Plus 2 5 °C—240 °C (260 °C); 30 m x 320 µm x 0.6 µm ---> Front Detector FID
	Aux EPC 1 He
	Aux EPC 2 He
	Aux EPC 3 He

Columns ▼ Options

Control Mode

	Actual	Setpoint
<input checked="" type="checkbox"/> On		
Flow	4 mL/min	<input type="text" value="4 mL/min"/>
Pressure	15.748 psi	<input type="text" value="15.749 psi"/>
Average Velocity	54.889 cm/sec	<input type="text" value="54.889 cm/sec"/>
Holdup Time	0.91093 min	<input type="text" value="0.91093 min"/>

Constant Flow ▼

Post Run:

Column #2 Configuration

((Initial): 0 min
 He @ 40 °C Oven
 Out: Ambient Pressure
 30 m x 320 µm x 0.6 µm

7.4.3 Oven

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns**
- Oven**
- Detectors
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous
 - Columns

Actual

Oven Temp On

40 °C

Equilibration Time:

Maximum Oven Temperature:

Override Column Max: 260 °C

	Rate °C/min	Value °C	Hold Time min	Run Time min
▶ (Initial)		40	3	3
*				

Post Run:

Post Run Time:

7.4.4 Detector

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors**
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous

FID

	Actual	Setpoint
<input checked="" type="checkbox"/> Heater:	250 °C	<input type="text" value="250 °C"/>
<input checked="" type="checkbox"/> Air Flow:	450 mL/min	<input type="text" value="450 mL/min"/>
<input checked="" type="checkbox"/> H2 Fuel Flow:	40 mL/min	<input type="text" value="40 mL/min"/>
<input checked="" type="checkbox"/> Makeup Flow: (N2)	40 mL/min	<input type="text" value="40 mL/min"/>
▶ Carrier Gas Flow Correction (None)		
<input checked="" type="checkbox"/> Flame	59.7 pA	
▶ No Column Comp		

7.4.5 Aux Heater

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors
- Aux Heaters**

Aux Heaters

Thermal Aux 2 (MSD Transfer Line)

Actual

On

250 °C

7.4.6 Signals

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous
 - Columns
 - Modules
 - ALS
 - Backflush
 - Readiness
 - GC Calculators

Dual	Signal Source	Data Rate / Min Peak Width	Zero	Save
F	#1: Front Signal (FID)	5 Hz / 0.04 min	Hz ? <input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
B	#2: None	50 Hz / 0.004 min	Hz ? <input type="checkbox"/>	<input type="checkbox"/>
B	#3: None	50 Hz / 0.004 min	Hz ? <input type="checkbox"/>	<input type="checkbox"/>
B	#4: None	50 Hz / 0.004 min	Hz ? <input type="checkbox"/>	<input type="checkbox"/>

Hide Dual Injection Signal Assignments

Signal Event Table

Signal Source	Time, min	Signal Event

Delete Events

7.4.7 Column Configuration

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous
 - Columns

Flow Paths:
 Front Inlet---->#1---->MSD
 Back Inlet---->#2---->Front Detector

Options

Column Outlet Pressure: 0 psi

Column	Calibration Results	Inlet	Outlet	Heated By
1 Restek 18004: 2037 RTX-BAC Plus 1 5 °C—240 °C (260 °C): 30 m x 320 µm x 1.8 µm	Uncalibrated	Front Inlet	MSD	Oven
2 Restek 18006: 2038 RTX-BAC Plus 2 5 °C—240 °C (260 °C): 30 m x 320 µm x 0.6 µm	Uncalibrated	Back Inlet	Front Detector	Oven
3 No Column Installed	Uncalibrated	Unspecified	Other	Oven
4 No Column Installed	Uncalibrated	Unspecified	Other	Oven

Options

Column Outlet Pressure: 0 psi

7.4.8 Module Configuration

Select...

- ALS
- ▼ Inlets
 - SSL - Front
 - SSL - Back
- Columns
- Oven
- Detectors
- Aux Heaters
- Events
- Signals
- ▼ Configuration
 - Miscellaneous
 - Columns
 - Modules
 - ALS
 - Backflush
 - Readiness
 - GC Calculators

Front Inlet
SS Inlet He

Back Inlet
SS Inlet He

Front Detector
FID
Makeup N2

Set Lit Offset with GC Keyboard.

Aux EPC 1,2,3
Aux EPC 1 He

Aux EPC 1,2,3
Aux EPC 2 He

Aux EPC 1,2,3
Aux EPC 3 He

7.4.9 GC Readiness

Select...	GC Readiness
ALS	Only checked components will affect the GC readiness
▼ Inlets	<input checked="" type="checkbox"/> Oven
SSL - Front	<input type="checkbox"/> Front Inlet (SS Inlet)
SSL - Back	<input checked="" type="checkbox"/> Back Inlet (SS Inlet)
Columns	<input checked="" type="checkbox"/> Front Detector (FID)
Oven	<input type="checkbox"/> Aux EPC 1
Detectors	<input type="checkbox"/> Aux EPC 2
Aux Heaters	<input type="checkbox"/> Aux EPC 3
Events	<input checked="" type="checkbox"/> Thermal Aux 2 (MSD Transfer Line)
Signals	

8 DATA ANALYSIS

8.1 Decision Criteria

The following criteria are applied through automated data analysis via Agilent MassHunter software. Integration parameters may be adjusted to effect proper integration. See TOX-101 for further guidance.

8.1.1 FID Method

8.1.1.1 Integration Criteria

Analyte	RT	%RT	Criteria	Integrator	Peak Filter
Methanol	1.227	10	Close RT	Agile2	Area ≥ 3000 counts
Ethanol	1.477	10	Close RT	Agile2	Area ≥ 3000 counts
Acetone	1.593	10	Close RT	Agile2	Area ≥ 50000 counts
Isopropanol	1.683	10	Close RT	Agile2	Area ≥ 3000 counts
T-butanol	1.857	10	Close RT	Agile2	Area ≥ 3000 counts

8.1.1.2 Calibration Criteria

Analyte	Curve Fit	Origin	Weight	Accuracy (+/-)	Levels (g/100mL)
Methanol	Linear	Ignore	1/x ²	15%	0.010, 0.025, 0.050, 0.100, 0.200, 0.400
Ethanol	Linear	Ignore	1/x ²	10%	0.010, 0.025, 0.050, 0.100, 0.200, 0.400, 0.500
Acetone	Linear	Ignore	1/x ²	15%	0.010, 0.025, 0.050, 0.100, 0.200, 0.400
Isopropanol	Linear	Ignore	1/x ²	15%	0.010, 0.025, 0.050, 0.100, 0.200, 0.400

8.1.1.3 Control Criteria

Analyte	Accuracy (+/-)	Levels
Methanol	15%	Cliniqa 1, 2
Ethanol	10%	Cliniqa 1, 2
Acetone	15%	Cliniqa 1, 2
Isopropanol	15%	Cliniqa 1, 2

8.1.1.4 Reporting Criteria

Analyte	Calculated Limit of Detection	Limit of Quantitation	Reporting Limit (Administratively Set)
Methanol	0.0019	0.010	0.005
Ethanol	0.0021	0.010	0.005
Acetone	0.0008	0.010	0.005
Isopropanol	0.0016	0.010	0.005

8.1.2 MSD Method

8.1.2.1 Integration Criteria

Analyte	RT	%RT	Criteria	Integrator	Quantifier Ion (m/z)	Peak Filter
Methanol	1.592	2	Close RT with Qualifiers	Agile2	31	Area ≥ 3000 counts
Ethanol	2.010	2	Close RT with Qualifiers	Agile2	31	S/N ≥ 10
Isopropanol	2.456	2	Close RT with Qualifiers	Agile2	45	S/N ≥ 10
Acetone	2.668	2	Close RT with Qualifiers	Agile2	43	S/N ≥ 10
T-butanol	2.909	2	Close RT with Qualifiers	Agile2	59	S/N ≥ 10

8.1.2.2 Qualifier Ion Criteria

Analyte	Qualifier Ion (m/z)	Relative Response	Criteria (+/-)
Methanol	29	72.0	24%
	32	79.0	37%
Ethanol	46	31.0	15%
	45	76.0	15%
Isopropanol	43	25.0	28%
Acetone	42	7.7	15%
	58	39.9	15%
t-butanol	57	10.0	15%
	41	23.0	30%

8.1.2.3 Calibration Criteria (Semi-Quantitative)

Analyte	Curve Fit	Origin	Weight	Accuracy (+/-)	Levels (g/100mL)
Methanol	Linear	Ignore	1/x ²	10%	0.010, 0.200
Ethanol	Linear	Ignore	1/x ²	10%	0.010, 0.200
Acetone	Linear	Ignore	1/x ²	10%	0.010, 0.200
Isopropanol	Linear	Ignore	1/x ²	10%	0.010, 0.200

8.1.2.4 Control Criteria

Analyte	Accuracy (+/-)	Level (g/100mL)
Methanol	10%	0.010
Ethanol	10%	0.010
Acetone	10%	0.010
Isopropanol	10%	0.010

8.1.3 Batch Acceptance

8.1.3.1 *Control Criteria*

Target analyte(s) shall not be detected in the Negative Control. Positive Control(s) shall have all target analytes identified. The software will automatically flag any control values that fail to meet the conditions in Section 8.

8.1.3.2 *Internal Standard*

The internal standard shall be recovered for all samples. The software will automatically flag any samples that exceed 10% variation in response of the calculated mean of the calibrators for that batch.

8.2 Calculations

8.2.1 MSD Screening

Calibration is linear with $1/x^2$ weighting. A two point semi-quantitative curve provides an estimated analyte concentration. For additional guidance, refer to Section 8.1.2.3 and TOX-101.

8.2.2 FID Confirmation

Calibration is linear with $1/x^2$ weighting. A six or seven point calibration curve is used to provide quantitative results (analyte dependent). Case samples are analyzed in duplicate and the values are averaged. For additional guidance, refer to Sections 8.1.1.2 and TOX-101.

8.2.3 Characterization of Whole Blood Controls

For commercial volatiles controls, each newly acquired lot of control will be analyzed at least 20 times in a minimum of four batches. The initial target value for the new control will be the average of these 20 values. Should the initial calculated value of the control exceed $\pm 5\%$ of the nominal value for ethanol (or ± 0.005 g/100mL, whichever is greater), or $\pm 10\%$ of the nominal value for any of the other volatiles, troubleshooting will occur before the control will be used in casework.

At least every six months, the accepted target value will be recalculated as the average value from all runs to date, excluding any failed analytical runs. Should the initial or recalculated target value of the control ever exceed $\pm 5\%$ of the nominal value for ethanol (or ± 0.005 g/100mL, whichever is greater), or $\pm 10\%$ of the nominal value for any of the other volatiles, the control may be degrading and a new lot should be purchased and characterized. The Technical Leader will ensure that a database of the lot performance of each lot of volatiles control is maintained.

9 REPORTING

9.1 Measurement Uncertainty

Refer to CHEM-100 and TOX-101.

TOX-200-20: Common Volatiles	Page 20 of 23	Issue Date: 03/06/2023
------------------------------	---------------	------------------------

9.2 MSD Screening

Analytes that are identified above the estimated 0.010 g/100mL reporting limit are confirmed by FID quantitative analysis prior to reporting. If no analytes are identified, then the results are reported as not detected.

9.3 FID Confirmation

Analytes are reported according to the following scheme:

Scenario	
Quantitated \geq 0.010 g/100mL	[analyte]: [concentration] [expanded measurement uncertainty]
Quantitated \geq 0.005 < 0.010 g/100mL	[analyte]: less than 0.010 g/100mL
Quantitated <0.005 g/100mL	[analyte]: not detected
Quantitated > highest calibrator	[analyte]: > [highest calibrator] g/100mL

9.4 Quantitative Values

Replicate values are averaged. This average value is truncated to three decimal places. The method's expanded uncertainty value is rounded up to the third decimal place. The current *k* value and a truncated coverage probability of 99.7% are also expressed.

Example:

Ethanol identified: 0.051 +/- 0.006; Acetone identified: 0.097 +/- 0.012; Methanol: not detected; Isopropanol: not detected; reported units g/100mL (grams per 100 milliliters). Measurement uncertainty provided at the 99.7% confidence level, $k=3.041$

10 CORRECTIVE MEASURES

If any criteria listed in Section 8 are not met, some or all of the following action steps may be appropriate (refer to TOX-101 for additional potential responses to QC failure(s)):

- Not reporting results from the batch and/or affected case samples
- Reaccession and reanalysis of the batch and/or affected case samples
- Performing instrument maintenance
- Remaking or using new reagents, calibrators, or control materials
- Notifying the Technical Leader who will ensure the root cause is determined and appropriate actions taken to address the issue(s)

11 PERFORMANCE CHARACTERISTICS

11.1 LOD (MSD and FID)

Analyte	FID LOD Calculated (g/100mL)	MSD Reporting Limit (g/100mL)
Ethanol	0.0021	0.010
Methanol	0.0019	0.010
Acetone	0.0008	0.010
Isopropanol	0.0016	0.010

11.2 LOQ (FID)

Analyte	Calculated (g/100mL)	Quantitation Reporting Limit (g/100mL)
Ethanol	0.0065	0.0100
Methanol	0.0057	0.0100
Acetone	0.0025	0.0100
Isopropanol	0.0049	0.0100

11.3 Linearity (FID)

Analyte	(g/100mL)
Ethanol	0.010 – 0.500
Methanol	0.010 – 0.400
Acetone	0.010 – 0.400
Isopropanol	0.010 – 0.400

11.4 Precision (FID)

Analyte	Low (%)	High (%)
Ethanol	1.83	1.60
Methanol	1.92	1.81
Acetone	5.83	5.13
Isopropanol	2.15	1.70

11.5 Carryover

There was no significant carryover for the MSD or FID methods at 0.400 g/100mL (methanol, acetone, isopropanol) or 0.500 g/100mL (ethanol).

11.6 Processed Sample Stability

When secured in unanalyzed, sealed headspace vials, samples are stable for at least 15 days in refrigerated conditions. Once the septum on a vial is punctured, the analyte response will decrease, becoming less stable after 24 hours. Samples may be reanalyzed for up to 24 hours after the initial analysis for screening purposes (GC-MSD analysis). Quantitative analyses (GC-FID) will not be reanalyzed.

12 LIMITATIONS

12.1 Interferences

No interferences have been identified for this method.

12.2 Interpretation

Ethanol is normally present in the human body at low levels (<0.001 g/100mL) due to bacterial fermentation in the intestines. Ethanol can also be produced because of putrefactive processes,

attributed to post-mortem processes and/or sample storage conditions. Consequently, caution should be exercised in the interpretation of low ethanol results (<0.04 g/100mL) in post-mortem cases.

13 SAFETY

Take standard precautions for the handling of chemicals and biological materials. Refer to the *FBI Laboratory Safety Manual* for guidance.

14 REVISION HISTORY

Revision	Issued	Changes
18	06/15/2022	Document reformat. Additional performance characteristics described in Section 11 .
19	02/01/2023	Updated storage conditions 5.6.1. Updated when to add an additional positive control 6.1.3 and 6.2.3.
20	03/06/2023	Clarified language in 3. Added storage tubes to 5.2. Clarified calibrator usage in 5.6.1. Added control storage guidance in 5.6.2. Clarified batch naming convention in 6.1. Updated 6.1.1 to be consistent with 6.2.1. Clarified batch naming convention and added optional intermediary container to 6.2. Added additional guidance for homogenizing samples and removed the option to dilute samples in 6.3. Updated Module Configuration screen capture in 7.2.7. Updated Gerstel AutoSampler Settings (FID and MSD) screen capture in 7.3. Updated Signals screen capture in 7.4.6. Updated qualifier ion criteria in 8.1.2.2. Clarified language used in 8.2.3.