

# Direct Solvent Extraction of Acid/Neutral Drugs from Biological Fluids

## Table of Contents

<b>1</b>	<b>INTRODUCTION .....</b>	<b>3</b>
<b>2</b>	<b>SCOPE.....</b>	<b>3</b>
<b>3</b>	<b>PRINCIPLE .....</b>	<b>3</b>
<b>4</b>	<b>SPECIMEN CRITERIA .....</b>	<b>3</b>
<b>5</b>	<b>EQUIPMENT .....</b>	<b>3</b>
5.1	Equipment .....	3
5.2	Consumables.....	3
5.3	Instruments .....	4
5.4	Software .....	4
5.5	Chemicals/Reagents .....	4
5.5.1	<i>Purchased .....</i>	<i>4</i>
5.5.2	<i>Prepared .....</i>	<i>4</i>
5.6	Standards/Controls .....	5
5.6.1	<i>Purchased .....</i>	<i>5</i>
5.6.2	<i>Prepared .....</i>	<i>5</i>
<b>6</b>	<b>PROCEDURE .....</b>	<b>7</b>
<b>7</b>	<b>ANALYTICAL PARAMETERS .....</b>	<b>9</b>
7.1	Agilent Gas Chromatograph.....	9
7.1.1	<i>Oven – Standard Conditions .....</i>	<i>9</i>
7.1.2	<i>Oven – Volatiles Analysis Conditions.....</i>	<i>9</i>
7.1.3	<i>Inlet/Carrier/Column .....</i>	<i>9</i>
7.2	Agilent Mass Spectrometer.....	9
<b>8</b>	<b>DATA ANALYSIS .....</b>	<b>11</b>
8.1	Decision Criteria.....	11
8.1.1	<i>Batch Decision Criteria.....</i>	<i>11</i>
8.1.2	<i>Chromatography .....</i>	<i>11</i>
8.1.3	<i>Mass Spectrometry.....</i>	<i>11</i>
<b>9</b>	<b>REPORTING .....</b>	<b>11</b>
<b>10</b>	<b>CORRECTIVE MEASURES .....</b>	<b>12</b>
<b>11</b>	<b>PERFORMANCE CHARACTERISTICS .....</b>	<b>12</b>
11.1	LOD.....	12
11.2	Carryover .....	13
<b>12</b>	<b>LIMITATIONS .....</b>	<b>13</b>
12.1	Interferences.....	13

<b>13</b>	<b>SAFETY.....</b>	<b>13</b>
<b>14</b>	<b>REVISION HISTORY.....</b>	<b>13</b>

## Direct Solvent Extraction of Acid/Neutral Drugs from Biological Fluids

### 1 INTRODUCTION

This procedure detects common acidic and neutral drugs in biological fluids.

### 2 SCOPE

Analyses	<input checked="" type="checkbox"/> Screening <input checked="" type="checkbox"/> Confirmation <input type="checkbox"/> Quantitation
Matrices	Biological fluids, food, beverages
Analytes	Selected acid-neutral drugs
Personnel	This document applies to authorized personnel who perform the described tasks, singly or in combination.

### 3 PRINCIPLE

Specimens are mixed with an internal standard, adjusted to an acidic pH, and extracted with an ether:toluene mixed solvent. Following centrifugation, the organic solvent is taken to dryness and the residue is partitioned between ethanol and hexane. The ethanol layer is taken to dryness and the extract is reconstituted in a chloroform/methanol mixture prior to analysis by GC/MS.

### 4 SPECIMEN CRITERIA

This procedure uses a biological fluid such as: blood, serum, plasma, urine, vitreous humor, or a prepared tissue homogenate. When available, 0.5 mL of blood or other fluids are used. This procedure may also be used to screen food and beverage samples for acidic and neutral drugs, provided that appropriate controls are simultaneously analyzed. A 0.5 g sample of a food or beverage homogenate or dilution is suggested for analysis.

### 5 EQUIPMENT

#### 5.1 Equipment

- A. Vortex mixer
- B. Centrifuge
- C. Rotator
- D. Evaporator with nitrogen

#### 5.2 Consumables

- A. 16 x 125 mm screw-top tubes with Teflon insert caps
- B. 16 x 100 mm culture tubes with polypropylene snap-tops
- C. 10 x 75 mm and 12 x 75 mm culture tubes with polypropylene snap-tops
- D. Routine laboratory supplies, including disposable pipettes, wooden sticks, test tube racks, graduated cylinders, etc.
- E. 10 cc glass centrifuge tubes (with conical bottom)
- F. ALS (automatic liquid sampler) vials 12x32 mm

### 5.3 Instruments

- A. Gas Chromatograph / Mass Spectrometer equipped with a 30 m x 0.25 mm x 0.25  $\mu$ m Rtx-5MS (or equivalent) column

### 5.4 Software

Component	Software	Version
Operating System	Microsoft Windows	7 Pro SP 1
GC/MS	Enhanced Chemstation	F.01.03.2357
Autosampler	Gerstel Maestro 1	1.5.3.3/3.5
Data Analysis	Xcalibur	4.0

### 5.5 Chemicals/Reagents

#### 5.5.1 Purchased

- A. N- Hexane (95% or equivalent)
- B. Potassium phosphate, monobasic (ACS grade or equivalent,  $\text{KH}_2\text{PO}_4$ )
- C. Diethyl ether (High purity grade or equivalent)
- D. Toluene (HPLC grade or equivalent)
- E. Chloroform ( $\text{GC}^2$  grade or equivalent)
- F. Methanol (Optima,  $\text{GC}^2$  grade or equivalent)
- G. Ethanol (Pharmaceutical grade or equivalent)
- H. Water (Deionized)
- I. Dichloromethane (Optima grade or equivalent)

#### 5.5.2 Prepared

**A. Potassium Phosphate Buffer Monobasic (5% w:v, pH 4.5):**

To a 100-mL volumetric flask, add 80 mL deionized water. Add 5 g [monobasic potassium phosphate](#) and mix well to dissolve. Bring to volume with deionized water, and verify  $4.0 < \text{pH} < 5.0$ . Store refrigerated in glass. Stable 1 month.

**B. Ether:Toluene (1:1 v:v):**

Combine 50 mL HPLC grade toluene with 50 mL diethyl ether. Mix well. Store in glass at room temperature. Stable 1 month.

**C. Chloroform:Methanol ( $\text{CHCl}_3$ :MeOH ) (4:1 v:v):**

Combine 40 mL chloroform with 10 mL methanol. Mix well. Store in brown glass at room temperature. Stable 1 month.

**D. Ethanol 80% (v/v aqueous):**

Measure 80 mL pharmaceutical grade ethanol into a 100-mL graduated cylinder. Bring to volume with deionized water and mix well. Store in glass at room temperature. Stable for 6 months.

## 5.6 Standards/Controls

### 5.6.1 Purchased

Storage and stability determined by manufacturer, unless otherwise noted.

**A. Methylphenylhydantoin (MPH)**

**B. Octanoic Acid**

**C. Negative Control:**

Purchased from Cliniq or an equivalent supplier, or prepared in-house from an appropriate blank specimen. Store refrigerated or obtain fresh.

**D. Barbiturate Mix-5:**

A mixture of five barbiturates at 250 µg/mL in methanol. Purchased from Cerilliant or another approved supplier. Contains amobarbital, butalbital, pentobarbital, phenobarbital, and secobarbital. This mixture may also be prepared from individual analyte stock solutions if necessary.

**E. Positive Control Solution Components:**

In addition to the Barbiturate Mix-5, target analytes (carbamazepine, carisoprodol, ibuprofen, meprobamate, and phenytoin) are obtained from an approved vendor in liquid (1 mg/mL) or solid form.

**F. Acetaminophen Standard (1 mg/mL):**

Purchased as a 1 mg/mL solution in methanol from Cerilliant or another approved supplier.

**G. Valproic Acid Standard (1 mg/mL):**

Purchased as a 1 mg/mL solution in methanol from Cerilliant or another approved supplier.

### 5.6.2 Prepared

**A. Methylphenylhydantoin Stock Standard (1 mg/mL):**

Add 10.0 mg of methylphenylhydantoin to a 10-mL volumetric flask. Dilute to the mark with methanol and mix well. Store refrigerated in glass. Stable for at least 2 years.

**B. Methylphenylhydantoin Working Internal Standard (30 µg/mL):**

Dilute 0.75 mL of the MPH Stock Standard to 25 mL with deionized water. Store refrigerated in glass. Stable for at least 2 years.

**C. Octanoic Acid Working Internal Standard (1mg/mL):**

Add 10.0 mg of octanoic acid to a 10-mL volumetric flask. Dilute to the mark with methanol and mix well. Store refrigerated in glass. Stable for at least 2 years.

**D. Positive Control Solution**

Solid analytes are dissolved in methanol or another appropriate solvent to prepare 1 mg/mL stock solutions. Analyte stock solutions (1 mg/mL) are added to a 25 mL volumetric flask which is brought to the mark with methanol as described in the table below:

Analyte(s)	Stock Conc. (mg/mL)	Spike Aliquot (μL)	Solution Volume (mL)	Solution Conc. (μg/mL)	Control Spike Aliquot (μL)	Matrix Volume (mL)	Control Conc. (ng/mL)
Barbiturate Mix-5	0.25	250	25	2.5	100	0.5	500
Carbamazepine	1	63	25	2.52	100	0.5	504
Carisoprodol	1	63	25	2.52	100	0.5	504
Ibuprofen	1	300	25	12	100	0.5	2400
Meprobamate	1	63	25	2.52	100	0.5	504
Phenytoin	1	63	25	2.52	100	0.5	504

The Positive Control Solution is stored refrigerated in glass or plastic. Stable for at least two years. (Note: Other drugs or metabolites may be added to this mixture as dictated by case needs with sufficient validation and/or analysis of concurrent controls.)

#### E. Positive Control:

100 μL of the Positive Control Solution is added to 0.50 mL of the Negative Control matrix on the day of analysis. Optional: 50 μL of a 1 mg/mL acetaminophen standard can be added directly to the Positive Control as well (yields a 100 μg/mL concentration). Other positive controls preparations may be used as is appropriate.

#### F. Valproic Acid Working Solution (100 μg/mL)





Dilute 500 μL of the Valproic Acid Standard (1 mg/mL) to 5 mL in methanol.

#### G. Volatiles Positive Control (20 μg/mL):

100 μL of the Valproic Acid Working Solution (100 μg/mL) is added to 0.50 mL of the Negative Control.

## 6 PROCEDURE

Step	Note	Reference/Lot
<b>A. Samples</b>		
1. To labeled 16 x 125 mm screw-top tubes add:		
<input type="checkbox"/> i. 0.5 mL of biological fluid		
<input type="checkbox"/> ii. 0.5 g of a prepared food homogenate		
<input type="checkbox"/> iii. 0.5 mL of a prepared beverage dilution		
<input type="checkbox"/> iv. 1 g of prepared tissue homogenate		
<b>B. Controls</b>		
<input type="checkbox"/> 1. Prepare Negative Control(s)		
<input type="checkbox"/> 2. Prepare Positive Control(s)		
<input type="checkbox"/> i. Prepare a Volatiles Positive Control if performing volatiles analysis		
<input type="checkbox"/> <b>C. QS to 1 mL with deionized water</b>		
<b>D. Internal Standard(s)</b>		
<input type="checkbox"/> 1. Add 25 µL of MPH Internal Standard Solution		
<input type="checkbox"/> i. Tissue specimens: add 0.5 mL of MPH Internal Standard Solution		
<input type="checkbox"/> ii. Volatiles analysis: add only 25 µL of Octanoic Acid Working Internal Standard Solution (1 mg/mL)		
<b>E. Buffer</b>		
<input type="checkbox"/> 1. Add 1 mL of 5% <a href="#">KH<sub>2</sub>PO<sub>4</sub></a> buffer solution		
<input type="checkbox"/> 2. Check pH to ensure pH is between 4 and 6		
<b>F. Extract</b>		
<input type="checkbox"/> 1. Add 5 mL of ether:toluene (1:1)		
<input type="checkbox"/> 2. Extract for 20 minutes on a rotator		
<input type="checkbox"/> 3. Centrifuge 5 minutes		
i. If emulsions develop, break up with wooden stick and recentrifuge		
<input type="checkbox"/> 4. Transfer organic (top) layer to a 16 x 100 mm tube		
<b>G. Concentrate/Clean Up</b>		
<input type="checkbox"/> 1. Evaporate solvent to dryness under nitrogen at 50°C.		

<input type="checkbox"/>	i. Volatiles Analysis: Pause evaporation at ~0.5 mL and skip to step I.		
<input type="checkbox"/>	2. Transfer dried residue to a 10 mL glass centrifuge tube with two successive 1 mL washes of CHCl <sub>3</sub> :MeOH (4:1)		
<input type="checkbox"/>	3. Evaporate solvent to a dry residue under nitrogen at 50°C.		
<input type="checkbox"/>	4. Add 2 mL hexane		
<input type="checkbox"/>	5. Add 200 µL 80% ethanol		
<input type="checkbox"/>	6. Vortex for 30 seconds		
<input type="checkbox"/>	7. Centrifuge for 5 minutes		
<input type="checkbox"/>	8. Discard the hexane (top) layer (optional)		
<input type="checkbox"/>	9. Transfer ethanol layer to a fresh 12 x 75 mm tube		
<input type="checkbox"/>	10. Evaporate ethanol layer to a dry residue under nitrogen at 50°C		
	<b>H. Reconstitute</b>		
<input type="checkbox"/>	1. Add 25-50 µL of CHCl <sub>3</sub> :MeOH (4:1)		
<input type="checkbox"/>	2. Vortex		
<input type="checkbox"/>	3. Transfer to ALS vial		
<input type="checkbox"/>	<b>I. Analyze 2 µL by GC-MS</b>		
	1. Analyze methylene chloride blanks at the sequence start and after positive control to reduce carryover		



## 7 ANALYTICAL PARAMETERS

### 7.1 Agilent Gas Chromatograph

#### 7.1.1 Oven – Standard Conditions

Step	Temperature (°C)	Hold (min)	Ramp (°C/min)
1	45	1	25
2	150	2	15
3	280	14	

Total Run Time (min): 29.87

#### 7.1.2 Oven – Volatiles Analysis Conditions

Step	Temperature (°C)	Hold (min)	Ramp (°C/min)
1	45	1	25
2	150	2	30
3	280	14	

Total Run Time (min): 25.53

#### 7.1.3 Inlet/Carrier/Column

Inlet	Carrier		Column		
Temperature (°C)	220	Gas	ultrapure helium	Type	DB-5MS
Injection Mode	Split	Mode	constant flow	Length (m)	30
Split Flow	12	Flow (mL/min)	1.2	Internal Diameter (mm)	0.25
Split Ratio	10:1			Film Thickness (µm)	0.25

### 7.2 Agilent Mass Spectrometer

	Standard Conditions	Volatiles Analysis Conditions
Ionization Mode	Electron Impact (+)	
Scan Mode	Full Scan	
Scan Range (m/z)	35-500	35-200
Solvent Delay (min)	3.0	3.5

Temperatures (°C)	
Source	230
Quadrupole	150
Transfer Line	280

## 8 DATA ANALYSIS

The authorized individual shall review electronic data files for acid/neutral drug screens and record that review in the exam documentation.

### 8.1 Decision Criteria

#### 8.1.1 Batch Decision Criteria

No analytes of interest will be detected in the Negative Control. For this purpose, analytes of interest are defined as those analytes that will be reported for this batch. All analytes should be detected in the Positive Control.

#### 8.1.2 Chromatography

The peak of interest will show good chromatographic fidelity, with reasonable peak shape, width, and resolution. In order to be determined acceptable, a chromatographic peak in an unknown sample will compare favorably to a chromatographic peak of the same analyte in a known sample analyzed on the same system in the same or subsequent analytical runs. Additionally, the following two criteria should be met.

##### 8.1.2.1 *Retention Time*

The retention time of the peak will be within  $\pm 2\%$  of the retention time (relative or absolute, as appropriate) obtained from injection of a reference standard, an extracted Positive Control, or an appropriate deuterated analog.

##### 8.1.2.2 *Signal-to-Noise*

To justify the existence of a peak, its baseline signal to peak-to-peak noise ratio will exceed 3. Further, the baseline signal for the peak of interest will be at least 10 fold greater than that for any observed peak at similar retention time in a Negative Control or solvent blank injected just prior to the sample.

#### 8.1.3 Mass Spectrometry

When necessary, the mass spectrum of the analyte of interest is compared to a reference standard or an extracted Positive Control. See the Guidelines for Comparison of Mass Spectra (TOX-104) for further guidance.

## 9 REPORTING

See Quality Control for Toxicology Examinations (TOX-101) for guidance on estimating the amount of an analyte in a specimen. When analyzing CAP T-Series or FTC specimens, if all decision criteria for an analyte of interest are met, but the concentration of butalbital, carbamazepine, carisoprodol, meprobamate, phenobarbital, phenytoin, and/or secobarbital is estimated to be below 1  $\mu\text{g/mL}$  (or 5  $\mu\text{g/mL}$  for acetaminophen) in two independent analyses, the analyte will not be reported. Note: the second analysis may be a repeat of this procedure or via another validated procedure. A Positive Control at the Cut-off Level is recommended for the second analysis.

TOX-201-09: Acid-Neutrals by GC/MS	Page 11 of 13	Issue Date: 02/11/2022
------------------------------------	---------------	------------------------

## 10 CORRECTIVE MEASURES

Refer to Quality Control for Toxicology Examinations (TOX-101) for guidance on action steps in the event of a quality control failure.

## 11 PERFORMANCE CHARACTERISTICS

### 11.1 LOD

Detection limits for common acidic and neutral analytes are listed in the table below. Note: LODs were not evaluated below 100 ng/mL, so true LODs for this method may be lower than what is listed in the table.

Analyte	Blood (µg/mL)	Urine (µg/mL)
Acetaminophen	25	1
Amobarbital	0.1	0.1
Brompheniramine	0.5	>0.5
Bupropion	0.1	0.25
Butalbital	2.5	0.5
Carbamazepine	0.1	0.1
Carisoprodol	0.5	0.1
Citalopram	0.25	0.5
Clozapine	>0.5	0.25
Cyclobenzaprine	0.1	0.1
Diphenhydramine	0.1	0.1
Ibuprofen	5	1
Ketamine	0.1	0.1
Lamotrigine	2.5	1
Lidocaine	0.25	0.25
Meprobamate	0.5	0.25
Methadone	0.25	0.1
Mirtazapine	0.1	0.1
Naproxen	50	2.5
Pentobarbital	0.1	0.1
Phenobarbital	0.1	0.1
Phenytoin	0.25	0.1

Propoxyphene	0.1	0.1
Secobarbital	0.1	0.1
Theophylline	2.5	0.1

## 11.2 Carryover

Carryover may occur after samples containing higher amounts of analyte; reinjection with appropriate solvent blanks may be performed.

## 12 LIMITATIONS

### 12.1 Interferences

None known. Grossly decomposed or putrefied samples, as well as samples that have been embalmed, may affect detection limits.

## 13 SAFETY

Take standard precautions for the handling of chemicals and biological materials. Refer to the FBI Laboratory Safety Manual for guidance.

## 14 REVISION HISTORY

Revision	Issued	Changes
09	02/11/2022	Document reformat. <u>1</u> -Simplified statement <u>2</u> -Reformatted scope statement <u>5</u> -Reorganized equipment listing, updated acetaminophen control amount <u>6</u> -Reformat of procedure for readability, eliminate form, specify injection amount <u>7</u> -Reformat GC/MS information Other minor edits and reorganization