

Physical Developer

Table of Contents

1	INTRODUCTION/SCOPE	2
2	STANDARDS AND CONTROLS	2
3	LIMITATIONS	2
4	EQUIPMENT	2
5	PROCEDURE	2
5.1	Solution Preparation	2
5.1.1	Maleic Acid solution.....	2
5.1.2	Tween 20 Solutions.....	2
5.1.3	DGME Solutions	3
5.2	Application	4
5.2.1	Tween 20 Application	4
5.2.2	DGME Application	4
5.3	Storage	4
5.4	Shelf Life.....	4
6	SAFETY	5
7	REVISION HISTORY	5

Physical Developer

1 INTRODUCTION/SCOPE

Physical Developer is used by FBI Laboratory Friction Ridge Discipline personnel to develop latent prints on porous and semi-porous surfaces.

2 STANDARDS AND CONTROLS

See *Processing Overview* (FRD-300).

3 LIMITATIONS

- A. All metal items, such as staples and paper clips, must be removed from item(s) prior to Physical Developer processing.
- B. Metal tweezers cannot be used during processing.

4 EQUIPMENT

- Distilled water
- Citric Acid
- Maleic Acid
- Silver Nitrate
- Tween 20
- Decaethylene Glycol Monododecyl Ether (DGME)
- Ferrous Ammonium Sulfate
- n-Dodecylamine Acetate
- Ferric Nitrate

5 PROCEDURE

5.1 Solution Preparation

Personnel will prepare the solutions as follows. Alternative amounts may be prepared, provided the same ratio of chemicals mixed is retained.

5.1.1 Maleic Acid solution

- A. Combine:
 - Maleic Acid - 25 g
 - Distilled water - 1000 mL
- B. Stir until solid dissolves.

5.1.2 Tween 20 Solutions

5.1.2.1 Tween 20 Redox solution

In the order listed below, stir until each solid is dissolved before adding the next solid:

- Distilled water - 1000 mL
- Citric Acid - 20 g
- Ferric Nitrate - 30 g

- Ferrous Ammonium Sulfate - 80 g

5.1.2.2 Tween 20 Detergent solution

- Combine:
 - n-Dodecylamine Acetate – 3 g
 - Synperonic N or Tween 20 – 4 g
 - Distilled water - 1000 mL
- Stir until all chemicals dissolve.

5.1.2.3 Silver Nitrate solution

- Combine:
 - Silver Nitrate - 200 g
 - Distilled water - 1000 mL
- Stir until solid dissolves.

5.1.2.4 Tween 20 Physical Developer working solution

- While stirring, combine in the order listed:
 - Redox solution - 1000 mL
 - Detergent solution - 40 mL
 - Silver Nitrate solution - 50 mL
- Stir for at least three minutes.

5.1.3 DGME Solutions

5.1.3.1 DGME Redox solution

Stir until each solid is dissolved before adding the next solid:

- Distilled water - 900 mL
- Citric Acid - 20 g
- Ferric Nitrate - 30 g
- Ferrous Ammonium Sulfate - 80 g

5.1.3.2 DGME Detergent solution

- Combine:
 - n-Dodecylamine Acetate – 1.5 g
 - DGME – 1.25 g
 - Distilled water - 1000 mL
- Stir until solid dissolves.

5.1.3.3 Silver Nitrate solution

- Combine:
 - Silver Nitrate - 200 g
 - Distilled water - 1000 mL
- Stir until solid dissolves.

5.1.3.4 DGME Physical Developer working solution

- While stirring, combine in the order listed:

- Redox solution - 900 mL
 - Detergent solution - 50 mL
 - Silver Nitrate solution - 50 mL
- B. Stir for at least ten minutes.

5.2 Application

5.2.1 Tween 20 Application

Personnel will complete the following steps in order:

- A. Immerse item(s) in Maleic Acid solution.
- B. Agitate solution, manually or with orbital shaker, for a minimum of 15 minutes.
- C. Immerse item(s) in Physical Developer working solution.
 - 1. Agitate solution, manually or with orbital shaker, for 10-15 minutes.
- D. Immerse item(s) in first water rinse for at least 1 minute.
- E. Rinse item(s) in second water rinse.
- F. Dry item(s) in air or by applying heat with an iron, heater, or dryer.
- G. Preserve appropriate friction ridge detail as applicable (e.g., digitally).
- H. Record the specific Physical Developer formula used in the case notes.

5.2.2 DGME Application

Personnel will complete the following steps in order:

- A. Immerse item(s) in Maleic Acid solution.
- B. Agitate solution, manually or with orbital shaker, for a minimum of 10 minutes or until bubbles no longer form in the solution (whichever takes longer).
- C. Immerse item(s) in Physical Developer working solution.
- D. Agitate solution, manually or with orbital shaker, for approximately 15 minutes.
 - 1. If latent print over-development begins to occur, samples may be removed prior to the 15 minutes.
- E. Immerse item(s) in first water rinse for at least 1 minute.
- F. Rinse item(s) in second water rinse.
- G. Dry item(s) in air or by applying heat with an iron, heater, or dryer.
- H. Preserve appropriate friction ridge detail as applicable (e.g., digitally).
- I. Record the specific Physical Developer formula used in the case notes.

5.3 Storage

- A. Maleic Acid, Redox, and Detergent solutions may be stored in any type of laboratory acceptable receptacle.
- B. Silver Nitrate solution must be stored in a dark bottle.
- C. Physical Developer working solution is not stored and is prepared as needed.

5.4 Shelf Life

- A. Maleic Acid and Redox solutions have indefinite shelf lives provided the reagent checks are satisfactory.

- B. Detergent and Silver Nitrate solutions have shelf lives of 1 year provided the reagent checks are satisfactory.
- C. Physical Developer working solution is not retained and is prepared as needed.

6 SAFETY

See FBI Laboratory Safety Manual for appropriate information.

7 REVISION HISTORY

Revision	Issued	Changes
04	07/01/2022	Format Updated. Section 1 – Scope updated.
05	10/15/2024	Synperonic N and Bleach rinse removed throughout and DGME solutions and application steps added. Distilled water added to equipment list. Requirement to record specific process used added.