# **Amido Black Fischer 98**

# **Table of Contents**

1	INTR	Introduction /Scope				
2	STANDARDS AND CONTROLS					
3	LIMI	2				
4	Equ	IPMENT	2			
5	Pro	CEDURE	2			
	5.1	Solution Preparation	2			
	5.1.	1 Standard Developer Solution	2			
	5.1.2 Alternate Developer Solution					
	5.2	Application				
	5.3 Storage					
	5.4	Shelf Life				
6	SAFE	ETY	4			
7	Revi	ISION HISTORY	4			

## **Amido Black Fischer 98**

## 1 Introduction / Scope

- A. Amido Black (Fischer 98) is a process used by FBI Laboratory Friction Ridge Discipline personnel to develop latent prints and enhance visible prints that have been deposited in blood.
- B. The process can be used on all surfaces but is primarily used on non-porous items.

#### 2 STANDARDS AND CONTROLS

See Processing Overview (FRD-300).

#### 3 LIMITATIONS

- A. The background of porous items may become stained during the process and obscure information.
- B. On clear or light colored surfaces, the alternate solution with Tween 20 may stain the background more than the standard version with Photo-Flo 600.

#### 4 EQUIPMENT

- Distilled water
- Water (for rinse)
- Naphthol Blue Black (dye content ≥85%)
- 5-Sulfosalicylic Acid (purity ≥99%)
- Formic Acid (concentrated)
- Sodium Carbonate
- Photo-Flo 600 Solution
- Tween 20
- N-dodecylamine Acetate
- Glacial Acetic Acid

#### 5 PROCEDURE

## 5.1 Solution Preparation

Personnel will prepare the solutions as follows. Alternative amounts of the final working solution may be prepared, provided the same ratio of chemicals mixed is retained.

# 5.1.1 <u>Standard Developer Solution</u>

- A. Combine the following:
  - o Naphthol Blue Black 3 g
  - o Glacial Acetic Acid 50 mL
  - o Distilled water 500 mL
  - o 5-Sulfosalicylic Acid 20 g
  - o Sodium Carbonate 3 g
  - o Formic Acid 50 mL
  - o Photo-Flo 600 Solution 12.5 mL

FRD-361-04: Amido Black Fischer 98	Page 2 of 4	Issue Date: 08/17/2022

- B. Stir solution until Naphthol Blue Black dissolves (approximately 30 minutes).
- C. Raise final volume to approximately 1000 mL with distilled water.
- D. Solution can be used immediately with acceptable results but works best if mixed and stored in a bottle several days before use.

## 5.1.2 <u>Alternate Developer Solution</u>

Photo-Flo 600 may be replaced with Tween 20 detergent solution.

#### 5.1.2.1 Tween 20 Detergent Solution

Personnel will combine the following and stir until all chemicals dissolve.

- n-Dodecylamine Acetate 3 g
- Tween 20 4 g
- Distilled water 1000 mL

# 5.1.2.2 Alternate Developer Solution

- A. Combine the following:
  - o Naphthol Blue Black 3 g
  - o Glacial Acetic Acid 50 mL
  - o Distilled water 500 mL
  - o 5-Sulfosalicylic Acid 20 g
  - o Sodium Carbonate 3 g
  - o Formic Acid 50 mL
  - o Tween 20 Detergent Solution 125 mL
- B. Stir solution until Naphthol Blue Black dissolves (approximately 30 minutes).
- C. Raise final volume to approximately 1000 mL with distilled water.
- D. Solution can be used immediately with acceptable results but works best if mixed and stored in a bottle several days before use.

## 5.2 Application

- A. Personnel will complete the following steps in order:
  - 1. Apply developer solution to the item by spraying, submersion, painting, or squirting.
    - Application can also be accomplished by the tissue method which involves wetting a durable tissue material and applying the material directly to the surface or by applying through a durable tissue material onto the surface.
  - 2. Leave developer solution on the item for 3 to 5 minutes.
  - 3. Rinse with water.
  - 4. Allow the item to dry.
- B. The developer solution may be reapplied as needed by repeating steps 1 through 3 until no further development is seen.
  - 1. Personnel will be cautious of overdevelopment and destruction of background.
- C. Capture appropriate friction ridge details as applicable (digitally or photographically).

FRD-361-04: Amido Black Fischer 98	Page 3 of 4	Issue Date: 08/17/2022
------------------------------------	-------------	------------------------

# 5.3 Storage

Developer solutions and Tween 20 Detergent solution may be stored in any type of laboratory accepted receptacle.

# 5.4 Shelf Life

Developer solutions and Tween 20 Detergent solution have an indefinite shelf life provided the reagent checks are satisfactory.

# 6 SAFETY

See FBI Laboratory Safety Manual for appropriate information.

## 7 REVISION HISTORY

Revision	Issued	Changes
03	07/15/2021	Replace Latent Print Units with Friction Ridge Discipline. Minor wording changes. Streamline equipment list. Change tissue to durable tissue material. Re-organization and re-numbering of sections.  Section 1 - removed blood fixer reference and added surfaces.  Section 2 - added limitations.  Section 4.1 - divided concept into two separate sections, Section 4.1.1 and Section 4.1.2 and added option on alternate amounts.  Section 4.2 - added reapplication allowance.  Section 5 - added Preamble reference.
04	08/17/2022	Reformatted  Section 5.1.2.1 – Separated Tween 20 Solution  Section 5.1.2.2 – Added full list for Alternate Developer Solution  Section 5.2 – Reworded tissue method  Section 5.3 – Added Tween 20 Solution  Section 5.4 – Added Tween 20 Solution