

Amido Black Fischer 98

Table of Contents

1	INTRODUCTION /SCOPE	2
2	STANDARDS AND CONTROLS	2
3	LIMITATIONS	2
4	EQUIPMENT	2
5	PROCEDURE	2
5.1	Solution Preparation	2
5.1.1	Tween® 20 stock solution	2
5.1.2	Amido Black Fischer 98 working solution	2
5.2	Application	3
5.3	Storage	3
5.4	Shelf Life.....	3
6	SAFETY	3
7	REVISION HISTORY	3

Amido Black Fischer 98

1 INTRODUCTION /SCOPE

Amido Black (Fischer 98) is used by FBI Laboratory Friction Ridge Discipline personnel to develop latent prints and enhance visible prints that have been deposited in blood on non-porous items.

2 STANDARDS AND CONTROLS

See Processing Overview (FRD-300).

3 LIMITATIONS

The process cannot be used on semi-porous and porous surfaces.

4 EQUIPMENT

- Distilled water
- Water (for rinse)
- Naphthol Blue Black (dye content $\geq 85\%$)
- 5-Sulfosalicylic Acid (purity $\geq 99\%$)
- Formic Acid (concentrated)
- Sodium Carbonate
- Tween® 20
- N-dodecylamine Acetate
- Glacial Acetic Acid

5 PROCEDURE

5.1 Solution Preparation

Personnel will prepare the stock and working solutions as follows: [Alternative amounts of the solutions may be prepared, provided the same ratio of chemicals mixed is retained.]

5.1.1 Tween® 20 stock solution

- A. Combine:
 - n-Dodecylamine Acetate - 3 g
 - Tween® 20 - 4 g
 - Distilled water - 1000 mL
- B. Stir until all chemicals dissolve.

5.1.2 Amido Black Fischer 98 working solution

- A. Combine:
 - Naphthol Blue Black - 3 g
 - Glacial Acetic Acid - 50 mL
 - Distilled water - 500 mL
 - 5-Sulfosalicylic Acid - 20 g

- Sodium Carbonate - 3 g
 - Formic Acid - 50 mL
 - Tween® 20 stock solution - 125 mL
- B. Stir solution until Naphthol Blue Black dissolves (approximately 30 minutes).
- C. Raise final volume to approximately 1000 mL with distilled water.
- D. Solution can be used immediately with acceptable results but works best if mixed and stored in a bottle several days before use.

5.2 Application

- A. Personnel will complete the following steps in order:
1. Apply Amido Black Fischer 98 working solution to the item
 2. Leave working solution on the item for 3 to 5 minutes.
 3. Rinse with water.
 4. Allow the item to dry.
- B. Working solution may be reapplied as needed using steps 1 through 4 above.
1. Personnel must be cautious of overdevelopment and destruction of background.
- C. Capture appropriate friction ridge detail as applicable (digitally or photographically).

5.3 Storage

Amido Black Fischer 98 working solution and Tween® 20 stock solution may be stored in any type of laboratory accepted receptacle.

5.4 Shelf Life

Amido Black Fischer 98 working solution and Tween® 20 stock solution have an indefinite shelf life provided the reagent checks are satisfactory.

6 SAFETY

See FBI Laboratory Safety Manual for appropriate information.

7 REVISION HISTORY

Revision	Issued	Changes
04	08/17/2022	Reformatted Section 5.1.2.1 – Separated Tween® 20 Solution Section 5.1.2.2 – Added full list for Alternate Developer Solution Section 5.2 – Reworded tissue method Section 5.3 – Added Tween® 20 Solution Section 5.4 – Added Tween® 20 Solution
05	06/03/2024	References to Photo Flo 600 removed throughout. Section 1 – Scope updated for substrate. Section 3 – Limitations updated for substrate and background information removed. Section 5.1 – solutions described as stock and working and standard developer solution removed with only Tween® 20 stock solution.

Revision	Issued	Changes
		Section 5.2 – specific application options removed.