

# Diaminobenzidine

## Table of Contents

<b>1</b>	<b>INTRODUCTION/SCOPE</b> .....	<b>2</b>
<b>2</b>	<b>STANDARDS AND CONTROLS</b> .....	<b>2</b>
<b>3</b>	<b>LIMITATIONS</b> .....	<b>2</b>
<b>4</b>	<b>EQUIPMENT</b> .....	<b>2</b>
<b>5</b>	<b>PROCEDURE</b> .....	<b>2</b>
5.1	Solution Preparation .....	2
5.1.1	Solution A (Fixer solution) .....	2
5.1.2	Solution B (Buffer solution) .....	2
5.1.3	Solution C (Diaminobenzidine solution) .....	2
5.1.4	Developer solution .....	3
5.2	Application .....	3
5.3	Storage .....	3
5.4	Shelf Life.....	3
<b>6</b>	<b>SAFETY</b> .....	<b>4</b>
<b>7</b>	<b>REVISION HISTORY</b> .....	<b>4</b>

# Diaminobenzidine

## 1 INTRODUCTION/SCOPE

- A. Diaminobenzidine is used by FBI Laboratory Friction Ridge Discipline personnel to develop latent prints and enhance visible prints deposited in blood.
- B. The process can be used on all surfaces but is primarily used on porous items.

## 2 STANDARDS AND CONTROLS

See *Processing Overview* ([FRD-300](#)).

## 3 LIMITATIONS

None

## 4 EQUIPMENT

- Distilled water
- 3,3'-Diaminobenzidine Tetrahydrochloride (purity  $\geq$  97%)
- Hydrogen Peroxide (30% Solution)
- 1 M Phosphate Buffer Solution (pH 7.4)
- 5-Sulfosalicylic Acid (purity  $\geq$  99%)

## 5 PROCEDURE

### 5.1 Solution Preparation

Personnel will prepare the solutions as follows. Alternative amounts of the final working solution may be prepared, provided the same ratio of chemicals mixed is retained.

#### 5.1.1 Solution A (Fixer solution)

- A. Combine:
  - 5-Sulfosalicylic Acid - 20 g
  - Distilled water - 1000 mL
- B. Stir until 5-Sulfosalicylic Acid dissolves.

#### 5.1.2 Solution B (Buffer solution)

- A. Combine:
  - 1 M Phosphate Buffer Solution - 100 mL
  - Distilled water - 800 mL

#### 5.1.3 Solution C (Diaminobenzidine solution)

- A. Combine:
  - 3,3'-Diaminobenzidine Tetrahydrochloride - 1 g
  - Distilled water - 100 mL
- B. Stir until 3,3'-Diaminobenzidine Tetrahydrochloride dissolves.

#### 5.1.4 Developer solution

- A. Combine:
  - o Solution B (Buffer solution) - 180 mL
  - o Solution C (Diaminobenzidine solution) - 20 mL
  - o Hydrogen Peroxide (30% solution) - 1 mL
- B. Mix thoroughly.

### 5.2 Application

- A. Personnel will complete the following steps in order:
  1. Apply Solution A (Fixer solution) to item by spraying, submersion, squirting, or painting.
    - i. Application can also be accomplished by the tissue method which involves wetting a durable tissue material and applying the material directly to the surface or by applying through a durable tissue material onto the surface.
    - ii. Leave on item for 3 to 5 minutes.
  2. Apply distilled water to item by spraying, submersion, squirting, painting, or tissue method.
    - i. Leave on item for 30 to 60 seconds.
  3. Apply Developer solution to item by spraying, submersion, squirting, painting, or tissue method.
    - i. Leave on item to achieve maximum contrast and development but do not exceed 5 minutes.
  4. Apply distilled water to item by spraying, submersion, squirting, painting, or tissue method to stop development process.
  5. Allow item to dry.
- B. Capture appropriate friction ridge details as applicable (digitally or photographically).

### 5.3 Storage

- A. Solution A (Fixer solution) and Solution B (Buffer solution) may be stored in any type of laboratory accepted receptacle.
- B. Solution C (Diaminobenzidine solution) must be stored frozen in a container that can withstand extreme cold.
- C. Storage is not applicable to developer solution. It is prepared as needed.

### 5.4 Shelf Life

- A. Solution A (Fixer solution) and Solution B (Buffer solution) have indefinite shelf lives provided the reagent checks are satisfactory.
- B. Solution C (Diaminobenzidine solution) has a 6 month shelf life when frozen, provided the reagent checks are satisfactory.
- C. Developer solution is not retained and is prepared as needed.

## 6 SAFETY

See [FBI Laboratory Safety Manual](#) for appropriate information.

## 7 REVISION HISTORY

Revision	Issued	Changes
02	07/15/2021	Replace Latent Print Units with Friction Ridge Discipline. Minor wording changes. Term “specimen” changed to “item” throughout. Streamline equipment list. Re-organization and re-numbering of sections. Section 1 - added last sentence. Section 3.1 - separated out into Section 3.1.1, Section 3.1.2, Section 3.1.3, and Section 3.1.4 and added ratio allowances for solutions. Section 4 - added Preamble reference.
03	08/17/2022	Reformatted <a href="#">Section 5.2</a> – Added description of tissue method