

General Microscopy Techniques

1 Scope

Personnel in the following categories of testing: hairs, fibers and textiles, glass, geologically-derived materials, and anthropology may use general microscopy techniques during their examination of evidence. Due to the wide variety of requests and examinations, this procedure provides general guidelines for general microscopy techniques and may or may not apply to every microscope used in each category of testing. This document applies to personnel within the Trace Evidence Unit (TEU) and Scientific and Biometrics Analysis Unit – Trace.

2 Equipment/Materials/Reagents

- Centering wrenches to fit microscope of interest
- Chemical resistant gloves
- Dust covers for microscopes
- Inert dusting gas
- Laboratory coat
- Lens cleaner
- Lens paper
- Petrographic microscope, with minimum 4x objective, 7.5x eyepieces rotatable stage, lower condensing lens
- Phase contrast microscope with minimum 10x magnification
- Prepared slide of particulate material immersed in mounting media with cover slip
- Sable brush
- Stereobinocular microscope, with minimum magnification of 4 diameters
- Water
- White polystyrene foam (small grain type)

3 Standards and Controls

Not applicable.

4 Sampling

Not applicable.

5 Procedure

5.1 General Care and Maintenance of Microscopes

The optical components of the microscopes should be kept clean at all times.

5.1.1 Replace the dust cover on the microscope at the end of each day, when available.

5.1.2 Clean the external surfaces of objectives, eyepieces and condensers.

5.1.2.1 Gently blow dust away, using a can of inert dusting gas, or using a clean soft sable brush.

5.1.2.2 After blowing away dust, use lens paper, dampened with water or lens cleaner to remove fingerprints, grease, oil and dirt. Never use alcohol on lenses.

5.1.3 White polystyrene foam (small-grain type) is recommended by some microscope companies for removing residues of immersion oil, skin grease and solvents. Break off a small piece and press a projecting part of it against the dry lens, rotating it as co-axially as possible with the objective. Any adhering foam granules can then be removed by blowing them away or using an absolutely clean sable brush.

5.1.4 Special care is required when working with acids and other chemical reagents. Their contact with the objectives should be strictly avoided. Clean the objectives (and other contaminated areas) at once after any accident. Even when particles are under a cover glass, there is a continuous stream of vapor from the corrosive substances that will impair the optical quality of the objective lens. Do not subject the lenses to prolonged exposure to these vapors.

5.2 Microscope Alignment

Modified Köhler illumination is used in petrographic and phase contrast microscopes.

5.2.1 Fully open field and aperture diaphragms, if present.

5.2.2 Adjust interpupillary distance (IPD).

5.2.3 Place specimen on stage.

5.2.4 Using only one eyepiece focus on the preparation using fine and coarse adjust controls on microscope. Adjust other eyepiece so that it is parfocal with the first.

5.2.5 For a Microscope Containing an Eyepiece Graticule

5.2.5.1 Without looking into the eyepieces, turn the eyelenses fully counterclockwise until they stop.

5.2.5.2 While looking into the eyepiece containing the graticule, adjust the eyepiece containing the graticule until the graticule comes into focus.

5.2.5.3 Do not adjust the eyepiece containing the graticule again.

5.2.5.4 Using the eyepiece containing the graticule, focus on the test object.

5.2.5.4.1 Place a flat object beneath the objective.

5.2.5.4.2 Using the lowest magnification and the eyepiece containing the graticule, bring the object into focus.

5.2.5.4.3 Using the highest magnification, bring the image into precise focus.

5.2.5.5 Adjust the eyepiece which does not contain the graticule.

5.2.5.5.1 Using the lowest magnification and the eyepiece which does not contain the graticule, focus the eyelens.

5.2.5.5.2 Using the highest magnification, look at the object with both eyepieces, and bring the image into precise focus.

5.2.6 For a Microscope that Does Not Contain an Eyepiece Graticule

5.2.6.1 Adjust the right eyepiece until it is in the center of its focusing range. Do not adjust the right eyepiece again.

5.2.6.1.1 Focus on the test object.

5.2.6.1.1.1 Place a flat object beneath the objective.

5.2.6.1.1.2 Using the lowest magnification and the right eyepiece bring the object into focus.

5.2.6.1.1.3 Using the highest magnification, bring the image into precise focus.

5.2.6.2 Adjust the left eyepiece

5.2.6.2.1 Without looking into the eyepiece, turn the eyelens fully counterclockwise.

5.2.6.2.2 Using the lowest magnification and the left eyepiece, focus the eyelens.

5.2.6.2.3 Using the highest magnification, look at the object with both eyepieces, and bring the image into precise focus, if necessary.

5.2.7 For a Microscope that Includes a Substage Assembly

5.2.7.1 Close field diaphragm by approximately one-half (1/2).

5.2.7.2 Focus and center condenser.

5.2.7.2.1 To focus the condenser, adjust the condenser focus knob until the image of the field diaphragm is in sharp focus.

5.2.7.2.2 To center the condenser, adjust the condenser centering screws until the image of the field diaphragm is centered in the field of view.

5.2.7.3 Open field diaphragm until just out of view.

5.2.7.4 Remove eyepiece and open aperture diaphragm until just out of view. Replace eyepiece.

5.2.8 For a Microscope with Adjustable Objectives

5.2.8.1 Place a prepared slide containing small particles on the microscope stage.

5.2.8.2 Using the lowest power objective, focus on a single particle.

5.2.8.3 Move the slide so that the particle of interest is in the center of field of view.

5.2.8.4 Rotate the stage. The particle should remain in the center of rotation. If the particle moves away from the center of the rotation:

5.2.8.4.1 Rotate the stage until the particle is furthest from the center of the field of view.

5.2.8.4.2 Adjust the objective centering wrenches until the particle is half way between its original position and the center of the field of view.

5.2.8.4.3 Repeat steps 5.2.8.4 through 5.2.8.4.2 until the particle remains in the center of the field of view upon rotation of the stage.

5.2.8.5 Repeat steps 5.2.8.2 through 5.2.8.4.2 for the remaining objectives.

5.2.9 For Microscopes with Phase Rings

5.2.9.1 If possible, align light source so that light from sub-stage is optimized while viewing through the microscope without a specimen on the stage.

5.2.9.2 Fully open the iris diaphragm.

- 5.2.9.3** With the two captive centering wrenches withdrawn, rotate the knurled turret ring to the open setting. These wrenches must be withdrawn to permit turret rotation.
- 5.2.9.4** Place a slide with mounted small particles on the stage.
- 5.2.9.5** Loosen condenser stop retaining screw and raise condenser with focusing knob.
- 5.2.9.6** Using the 10x objective, focus on the particles.
- 5.2.9.7** Center annular rings.
- 5.2.9.7.1** Place specimen slide on stage and open field diaphragm until iris just disappears from field of view.
- 5.2.9.7.2** Bring the specimen into good focus.
- 5.2.9.7.3** Turn the knurled ring to the 10x annulus and push in the two centering wrenches.
- 5.2.9.7.4** Image the aperture and bring the annulus of the condenser and the diffraction plate of the objective into sharp focus by turning the focusing knob of the aperture viewing unit.
- 5.2.9.7.5** Adjust the annular diaphragm centering wrenches until the annular image is superimposed on the diffraction plate.
- 5.2.9.7.6** Bring the specimen into sharp focus.

5.3 Color Balancing – Comparison Microscopes

- 5.3.1** After performing modified Köhler illumination on both sides of the comparison microscope, place slides with colored fibers from the same known sample on both sides of the comparison microscope.
- 5.3.2.** Adjust the light voltage (if applicable), aperture diaphragm, and field diaphragm until the background color appears the same on both sides of the comparison microscope.
- 5.3.3** If so equipped, adjust the Vari-Lum or the colored filters until the colors appear the same in both sides of the comparison microscope.

5.4 Microscope Transport

If items will be examined outside the Laboratory there may be a need to bring a stereobinocular microscope to the external location to facilitate the examinations.

5.4.1 The microscope will be disassembled into its component parts and carefully packed in a container (*e.g.*, Pelican case) to ensure they will not be damaged during transport.

5.4.2 The component parts will be inspected for possible damage and the microscope will be reassembled at the external location. If any damage has occurred, the microscope will not be used.

6 Calculations

Not applicable.

7 Measurement Uncertainty

Not applicable.

8 Limitations

Not applicable.

9 Safety

No potential hazards are associated with these techniques.

10 References

- Nikon, How To Use A Microscope And Take A Photomicrograph, Nikon Corporation, 1998.
- Patzelt, Walter J., Polarized Light Microscopy, 3rd edition, Ernst Leitz Wetzlar GmbH., 1985.
- Möllring, F.K. Microscopy From The Very Beginning, Carl Zeiss, Oberkochen, West Germany.
- Handbook of Incident Light Microscopy, Carl Zeiss, Oberkochen, West Germany.

- Determann, H. and F. Lepusch, The Microscope and its Application, Ernst Leitz Wetzlar GmbH.
- Delly, J.G., The Michel-Lévy Interference Color Chart-Microscopy's Color Key, *Microscope* 37, 89-102, 1989.

Rev. #	Issue Date	History
2	02/07/2018	Updated title removing reference to TEU to reflect change from unit-specific to discipline-specific. Added more background information and Scientific Analysis Unit - Trace to Scope in Section 1. Updated Sections 2, 5.1, 5.1.1, and 5.1.2.1 for clarity, Removed Section 4 Calibration. Updated new Section 4 by adding Sample Selection.
3	02/10/2020	Updated 'geological materials' to 'geologically-derived' in Scope. Updated Scientific and Biometrics Analysis Unit – Trace name. Removed 'Sample Selection' from Section 4 title. Updated language in Sections 5.3.1 and 5.3.2 for clarity. Added Section 5.4.

Approval

Redacted - Signatures on File

Trace Evidence Unit
 Chief

Date: 02/07/2020

Scientific and Biometrics
 Analysis Unit Chief

Date: 02/07/2020

Mineralogy Technical
 Leader

Date: 02/07/2020

Hair and Fiber Technical
 Leader

Date: 02/07/2020

Anthropology Technical
 Leader

Date: 02/07/2020

Trace Evidence Unit Evidence Processing Procedures

1 Scope

1.1 This document describes procedures used by personnel within the Trace Evidence Unit for the processing of physical evidence for trace evidence. This includes the hair, fibers and textiles, geologically-derived materials, glass, and anthropology categories of testing.

1.2 The nature and extent of processing will be determined by the type of evidence submitted, by previous handling or examinations of the evidence, and by the requested examinations of the contributor.

1.3 Guidelines for the microscopic analyses and comparison of hairs, textile fibers, anthropological and geologically-derived materials can be found in the category of testing's specific protocols.

2 Equipment/Materials/Reagents

- Stereobinocular microscope, magnification range from approximately 2x to 40x
- Permount mounting medium
- Xylene substitute, Xyless or Xylene
- Cargille liquids or equivalent
- Glass microscope slides and coverslips
- Kraft paper
- Pillboxes
- Forceps
- Spatula
- Scissors
- Probes
- un-du[®]
- Disposable lab outerwear
- Cavicide or equivalent cleaning solution
- Lint free wipe
- Paper, various sizes
- Rock hammer
- Mortar/pestle
- KPac or equivalent
- Glassine envelopes
- Lux-o-lamp

- Alternate Light Source
- Adjustable rack
- Blotter paper
- Single-use vacuum filters
- Vacuum
- Laboratory coats
- Gloves

3 Standards and Controls

Not applicable.

4 Sampling

Due to the condition of the submitted items and the number of hairs and/or fibers present, a single sample selection scheme cannot account for all possible scenarios. As a result, Examiner and/or Physical Scientist discretion will determine a representative sample on a case-by-case basis when all hairs and/or fibers are not mounted. A representative sample is a selection of hairs and/or fibers that captures the varying characteristics of the total hairs and/or fibers collected from an item of evidence.

For a known fiber sample selection, the representative sample will represent the range of colors and fiber types comprising the item.

5 Procedures

When items of evidence that are received in the TEU Hair and Fiber Group and listed on the assigned case record will not require examinations, the Examiner will document this decision in the Case Record Communication Log.

For cases containing more than one item of evidence, prior to processing any items of evidence that require the use of a processing room(s), a processing plan will be documented in the Case Record Communication Log by the TEU personnel who will be processing the evidence. This processing plan will be reviewed by any qualified TEU personnel prior to the start of processing. The review will be documented in the Case Record Communication Log.

5.1 Processing Physical Evidence

5.1.1 Before evidence is processed, the processing area and all utensils (*e.g.*, forceps and scissors) will be cleaned using at a minimum, a cleaner such as Cavicide and a lint free wipe.

5.1.2 Gloves, at a minimum, will be changed between cases. Other PPE will be changed as necessary. Facemasks will be worn during the processing of items/cases that have potential for DNA analysis.

5.1.3 All evidence will be processed over clean paper that is placed on the surface of the table.

5.1.3.1 A clean sheet of paper will be used for the processing of each item of physical evidence unless case circumstances indicate otherwise. Items received in the same packaging may be processed on the same piece of paper.

5.1.4 Accessory lighting, special lighting techniques, and magnification may be used as needed.

5.1.5 The item of evidence will be described regarding type, color, size, and style, and carefully evaluated to determine its condition including damage, stains, etc. The item will be marked with the item number and initials of the processor when possible.

5.1.6 Additional items recovered during processing will be documented in the notes and may be subdivided; *e.g.*, Item 1 Pants becomes:

- Item 1 Pants
- Item 1-1 Belt from Item 1

5.1.7 Items of evidence that are scraped will be either hung on an adjustable rack above a table or manually handled, depending on the size of the item. Smaller items such as knives, sticks, gloves, etc. may be processed at a workstation using a lux-o-lamp or stereobinocular microscope. When processing items at a workstation or stereobinocular microscope, work surfaces will be cleaned using at a minimum a cleaner (*e.g.*, Cavicide, cleaning wipe). Any tools used for processing evidence will be cleaned prior to processing.

5.1.8 Visible debris can be picked off of the item and preserved in a separate pillbox. Soil that may be layered will be collected to preserve the layer structure. Layered soil can be either picked off an item and placed in a separate pillbox or a cutting of the item containing the soil can be collected and placed in a separate pillbox.

5.1.9 The adjustable rack to which the item is attached will be adjusted to allow the item to hang just above the surface of the table.

5.1.10 The item will be gently scraped to remove trace evidence that is adhering to the surface of the item. As the case warrants, debris may also be recovered by vacuuming the object and collecting the debris on a vacuum filter, or by taking tape and patting it across the surface of an item.

5.1.11 Debris removed from the inside of items may be separated from outside debris, as warranted by the circumstances of the case.

5.1.12 When processing shoes for glass, after scraping the items, the soles of the shoes will be examined for the presence of embedded glass.

5.1.12.1 The soles of the shoes will be assessed for cuts or tear.

5.1.12.2 Insert a metal probe into the cut or tear. If any solid objects are embedded in the cut or tear, gently pry the object out, and place it in a separate pillbox.

5.1.13 Items recovered from pockets of submitted clothing may be placed in a separate packaging and appropriately marked with laboratory number, item number and the processor's initials.

5.1.14 Debris removed from an item may be either collected in a pillbox or other suitable container, or directly mounted on a glass microscope slide following the procedures outlined below. The receptacle will be appropriately marked with the Laboratory number, item number and initials of the processor. Pillboxes generated from evidence from different locations (e.g., Victim, Subject or Crime Scene) will be placed in separate bags.

5.1.15 If necessary, a known sample of fabric will be removed and placed in a druggist fold, appropriately marked with the Laboratory number, item number and initials of the processor. If necessary, the location of the sample site will be documented via diagrams, descriptions or other equivalent means in the case notes.

5.1.16 When the item of evidence has been processed, it will be returned to its original container and sealed. If the original container is replaced or damaged, a new one will be furnished, indicated in case notes and the original packaging will be retained within new packaging if possible.

5.1.17 After all items have been processed, they will be properly stored in a secure evidence cabinet, refrigerator, safe or bulky evidence room.

5.1.18 All evidence packages and/or boxes stored in any cabinet, refrigerator, safe or bulky storage that is not under active examination must be under proper seal. The Laboratory number should be clearly visible.

5.2 Procedures to be Used When Evidence Processing Includes Scraping Items Submitted from a Combination of Suspect(s), Victim(s), and Crime Scene(s)

5.2.1 To protect against cross-transfer and contamination during scraping, items submitted from the victim(s) and items submitted from the subject(s) will be processed in a different room or on a different date. The processing room used will be documented in the case notes. If items

are processed at an individual's workspace directly, the area will be thoroughly cleaned between victim, suspect and crime scene items.

5.2.2 A clean laboratory coat and protective gloves will be worn prior to entering the processing room.

5.2.3 The processing room will be thoroughly cleaned according to steps described in TEU Evidence Processing Procedure, sections 5.1.1 and 5.1.2.

5.2.4 When the processing of items from the victim is complete, the technician will remove their laboratory coat and discard protective gloves.

5.2.5 A clean laboratory coat and new protective gloves will be worn prior to entering the second processing room.

5.2.6 The processing room will be thoroughly cleaned according to steps described in TEU Evidence Processing Procedure, sections 5.1.1 and 5.1.2.

5.2.7 Clean utensils will be used to remove debris from items and for taking known fabric samples.

5.2.8 Items submitted from the crime scene may or may not be processed in the same room(s) as the victim/subject, depending on the circumstances of the case. Clean utensils and examination paper will be used in either case.

5.2.9 No specific order is required in processing of victim, subject, and crime scene items.

5.3 Debris Screening and Slide Preparation - Hairs

5.3.1 The screening of the evidentiary item or recovered debris is facilitated by the use of lux-o-lamp magnifiers and/or stereobinocular microscopes. This debris will be examined for the presence of hair evidence.

5.3.2 When a large number of hairs are present in the debris, a representative sample of hairs of different ancestral characteristics, body area, length, color, texture and thickness may be mounted. The number of hairs mounted on glass microscope slides may be influenced by the types of hairs in the questioned debris and the circumstances of the case. The letters "R/S" are placed on the glass microscope slide and included in the case notes to indicate that a representative sample of hairs was mounted and that additional hairs are present on or in the item/pillboxes.

5.3.2.1 If a representative sample of hairs from the questioned item(s) has been mounted on glass microscope slides, a targeted search of the questioned item(s) or recovered debris will be conducted if the following criteria are met:

- No hair association was found.
- A known hair sample(s) has been submitted from an individual(s) unrelated to the identified source of the questioned item(s) (*e.g.*, known hairs from victim, questioned items from suspect).

5.3.3 The screening of the debris can be directed by a “target” search, *e.g.*, searching for certain types of hairs (apparent head or pubic hairs) or searching for hairs that are similar to hairs comprising a known hair sample, when available.

5.3.4 When mounting several hairs of different lengths on a single slide, the length of the longest hair will be recorded on the frosted end of the glass microscope slide.

5.3.5 Hairs will be mounted on a clean glass microscope slide using a suitable mounting medium, *e.g.*, Permount. Each slide will contain the Laboratory number, the item number, and the initials of the individual preparing the slide.

5.3.6 Placing a thin film of solvent (such as Xylene substitute) on the surface of the slide will allow hairs to adhere temporarily until the mounting medium is applied. Using clean forceps, hairs will be placed on the slide and arranged so they can be completely covered by the glass coverslip.

5.3.7 Excess solvent will be blotted off to avoid run-off of the excess solvent when the mounting medium is applied and to help arrange hairs on the glass microscope slide. The used blotter paper will be discarded in the appropriate receptacle (see FBI Laboratory Safety Manual) between slides.

5.3.8 Forceps will be carefully cleaned between different items/pillboxes.

5.4 Debris Screening and Slide Preparation - Fibers

5.4.1 The screening of evidentiary items or recovered debris is facilitated by the use of lux-o-lamp magnifiers and stereobinocular microscopes. This debris will be examined for the presence of fiber evidence.

5.4.2 When the number of fibers present in the debris is such that all of the fibers cannot reasonably be mounted on a glass microscope slide, a representative sample of fibers of different colors, shapes and sizes may be mounted. The number of glass microscope slides prepared during the initial screening is dependent on the number and types of fibers in the questioned debris and the circumstances of the case. The letters “R/S” are placed on the glass microscope slide and included in the case notes to indicate that a representative sample of fibers was mounted and that additional fibers are present in the item/pillbox.

5.4.3 The screening of the debris can also be directed by a “target” search, *i.e.*, looking for fibers that are similar to fibers comprising a known fiber sample, when available. The known fiber sample may include carpet samples and/or fabric samples either submitted separately or collected during the processing of clothing items.

5.4.4 Fibers will be mounted on a clean glass microscope slide using a suitable mounting medium, *e.g.*, Permount. Each slide will contain the Laboratory number, the item number, and the initials of the individual preparing the slide.

5.4.5 Placing a thin film of solvent (such as Xylene substitute) on the surface of the slide will allow fiber samples to adhere temporarily until the mounting medium is applied. Using clean forceps, fibers will be placed on the slide and arranged so they can be completely covered by the glass coverslip.

5.4.6 Excess solvent will be blotted off to avoid run-off of the excess solvent when the mounting medium is applied and to help arrange fibers on the glass microscope slide. The used blotter paper will be discarded in the appropriate receptacle (see FBI Laboratory Safety Manual) between slides.

5.4.7 Forceps will be carefully cleaned between different items/pillboxes.

5.4.8 When complete yarns are identified in the pillbox debris, they will be characterized (*e.g.*, diameter, twist, construction) before being separated and mounted on the glass microscope slide. Consideration will also be given to physically matching yarns to damaged fabric before mounting fiber samples from the yarn on a slide.

5.5 Selection and Preparation of Known Fiber Slides

5.5.1 A known sample will be selected that represents the range of colors and fiber types comprising the textile.

5.5.2 If possible, known yarn samples will not be taken from damaged areas because of potential future yarn/fabric matches.

5.5.3 Fiber samples from yarn types present in the fabric will be mounted. Warp yarns and fill yarns may be separately mounted. Sewing thread and button thread fiber samples may also be mounted.

5.5.4 In addition to the Laboratory number, item number and the initials of the processor, the letters “kn”, which indicates a known sample, will be written on the frosted end of the glass microscope slide.

5.6 Debris Screening and Sample Preparation – Geologically-Derived Materials

5.6.1 Items are air-dried, if needed.

5.6.2 Individual components are removed from the items for identification as necessary. These individual sub-samples may not be representative of the entire item. Sub-samples are chosen based on the need to identify a particular component by a specific technique, and by its availability or presence in an item.

5.6.2.2 Items may be mechanically broken to facilitate sub-sampling.

5.6.3 Soil removed from objects, *e.g.*, shoes, will be kept as coherent as possible.

5.7 Debris Screening and Sample Preparation – Glass

The screening of evidentiary items or recovered debris is facilitated by the use of lux-o-lamp magnifiers and stereobinocular microscopes. This debris will be examined for the presence of glass evidence. As needed, refer to the *Forensic Glass Examinations, Section 4; Refractive Index of Glass by GRIM, Section 4; and Elemental Analysis of Glass by Inductively Coupled Plasma – Optical Emission Spectrometry (ICP-OES), Section 4; Laboratory Annealing of Glass, Section 4.2* procedures for sample preparation and sampling guidance for each analysis technique.

5.8 Envelope/Letter Processing

5.8.1 When envelopes/letters containing stamps or labels are received in the TEU, an attempt shall be made to remove the stamps or labels using a suitable solvent such as un-du®.

5.8.2 If the stamp or label cannot be removed using a suitable solvent, no further attempt will be made to remove the stamp or label.

5.8.3 If the stamp or label can be removed using a suitable solvent, the underside shall be examined for the presence of trace evidence. If trace evidence is identified, it shall be handled as described previously.

5.8.3.1 Each stamp or label that is removed shall be, subdivided and identified; *e.g.*, Item 1 Envelope becomes:

- Item 1 Envelope
- Item 1-1 Stamp from Item 1
- Item 1-2 Label from Item 1

5.8.3.2 Each stamp or label that is removed shall be placed on a suitable clear plastic sheeting.

5.9 Adhesive Surfaces

Items of evidence may contain adhesive surfaces that require other units within the laboratory to examine them prior to removal. If adhesive surfaces are observed and any affected unit(s) is already on the examination plan to examine the item(s), the affected unit will be contacted. This contact is to determine how far the trace evidence examination process can proceed without interfering with the examinations of the affected unit(s) examinations. This communication will be recorded in the case communication log. If the affected unit(s) is not on the examination plan, then the examiner or scientist will follow the appropriate Laboratory Operations Manual practice and contact the appropriate individual(s) to determine if this examination(s) is necessary and requested. If the item(s) needs to be returned for additional trace evidence examinations, the examiner or scientist will mark the outer packaging of the item(s) to indicate that the item(s) needs to be returned for additional trace evidence examinations prior to the adhesive item(s) being removed.

6 Calculations

Not applicable.

7 Measurement Uncertainty

Not applicable.

8 Limitations

Not applicable.

9 Safety

9.1 While working with physical evidence, Laboratory personnel will wear appropriate protective attire (at a minimum, a laboratory coat and gloves).

9.2 Universal precautions will be followed.

9.3 No specific hazards are associated with the microscopic examination techniques performed.

9.4 Refer to the safety data sheet (SDS) for guidelines regarding the use of a specific chemical.

10 References

- FBI Laboratory Quality Assurance Manual.
- FBI Laboratory Operations Manual.
- FBI Laboratory Safety Manual.
- Trace Evidence Procedures Manual.

Rev. #	Issue Date	History
3	10/02/2017	<p>Section 1: Updated Scope 1.1 and 1.2. Section 4: Updated wording for Sampling. Section 5.1.1: Updated wording. Section 5.1.2: Added section to have gloves changed between cases. Masks added for DNA cases. Section 5.2: Added entire section. Section 5.4.2 and 5.4.3: Edited sections to allow for the choice of target searches or representative samples based on case circumstances. Renumbered remainder of document due to added section. Section 5.10: Added section covering items with adhesive surfaces.</p>
4	02/10/2020	<p>Added wording to Section 5 and removed reference to Legacy cases. Section 5.6.2.1 deleted. Updated wording to Sections 5.4.8, 5.8, and 5.9.3.2 Changed 'geological' to 'geologically-derived' throughout. Combined Sections 5.6 and 5.7 into one section for Geologically-Derived Materials and made more general to apply to all processing. Renumbered remainder of document due to removed section.</p>

Approval

Redacted - Signatures on File

Trace Evidence Unit
 Chief

Date: 02/07/2020

Mineralogy Technical
 Leader

Date: 02/07/2020

Hair and Fiber Technical
 Leader

Date: 02/07/2020

Anthropology Technical
 Leader

Date: 02/07/2020