

FDDU

Procedures for STR Amplification

1 Scope

These procedures apply to DNA personnel who perform direct amplification of nuclear DNA with the polymerase chain reaction (PCR) using the AmpF/STR® Identifiler® Direct Amplification Kit or the GlobalFiler® Express PCR Amplification Kit on DNA databasing samples received by the Federal DNA Database Unit (FDDU).

2 Equipment/Materials/Reagents

Equipment/Materials

- General Laboratory Supplies
- STACS (Sample Tracking and Control System) Software (STACS DNA Inc.), version 6.0 or higher
- Barcode printer with appropriately sized labels (2.0" x 0.5" or equivalent)
- Barcode Scanner, Hand-held (Symbol LS4000i, 4008i, LS4071 or equivalent)
- Plate Sealer, microplate (Agilent PlateLoc or equivalent) with Heat Sealing Foil
- Optical Adhesive Covers (Applied BioSystems or equivalent)
- Tecan Robotic Workstation (Master Mix Addition)
 - Tecan EVOware Software, version 2.0 or higher (Tecan)
- Optical Cover Compression Pads (Applied BioSystems or equivalent)
- Thermal Cycler (Applied BioSystems GeneAmp® PCR System 9700 or ProFlex™ PCR System)
- Plate centrifuge (Eppendorf Model 5804 or equivalent)
- Laboratory incubator (Isotemp or equivalent)

Reagents

- Liquinox™ Critical Cleaning Liquid Detergent (Alconox or equivalent)
- Water, reagent grade (VWR or equivalent)
- Bleach, 3% (household or equivalent)
- 9947A, 10 ng/μl (Promega DD1001, DD100A or equivalent)
- AmpF/STR® Identifiler® Direct PCR Amplification Kit (Applied BioSystems)
 - Contains AmpF/STR® Identifiler® Direct Master Mix (includes AmpliTaq Gold® enzyme) and AmpF/STR® Identifiler® Direct Primer Set
- GlobalFiler® Express PCR Amplification Kit (Applied BioSystems)
 - Contains GlobalFiler® Express Master Mix, Master Mix Additive, GlobalFiler® Express Primer Set, GlobalFiler® Express Allelic Ladder, and DNA Control 007

3 Standards and Controls

Two Combo controls (aka negative) and two Blood/Buccal Internal Standard (BIS) controls are included on each amplification plate.

The BIS is prepared as described in the DNA QA procedure for reagents (i.e., DNAQA 609) and is added to the plate during the sample punch procedure. These controls will be interpreted according to the criteria in the applicable FDDU procedure (i.e., FDDU 315).

4 Procedures

Refer to the DNA Procedure Introduction (i.e., DNAQA 600) for applicable laboratory quality assurance and cleaning instructions.

When using a Robotic Workstation, ensure general instrument cleaning and maintenance is done prior to use, as needed. See Appendix A for additional guidance.

4.1 Processing with SwabSolution

For plates prepared using SwabSolution, perform the following steps any time prior to loading the plate onto the Tecan instrument or before manually adding master mix to the plate.

- 4.1.1 Ensure that the plate is adequately sealed after punching, then quick-spin the 96-Well plate(s) in a centrifuge (approximately 30 seconds).
- 4.1.2 Incubate the plate(s) at 70°C for 10 minutes.
- 4.1.3 After incubation, quick-spin the 96-Well plate(s) again in a centrifuge (approximately 30 seconds) before proceeding to Master Mix addition.

4.2 Master Mix (MM) Addition

The Master Mix Addition procedure may be performed based on the amplification kit and sample type as follows:

	Identifiler Direct MM	GlobalFiler Express MM
Blood	N/A	Manual or Automated
Buccal	Manual or Automated	

- 4.2.1 Prior to the first use of a GFE kit, prepare the GFE Master Mix tube by adding Master Mix Additive to the GFE Master Mix tube.
 - For a 200 reaction kit: Add 80µL of Additive
 - For a 1000 reaction kit: Add 390µL of Additive

NOTE: Amplification WILL FAIL without the Master Mix Additive in the GFE Master Mix.

- 4.2.2** Prepare the amplification master mix for the required number of samples. Include extra samples in the calculation for overage. Ensure the preparations of the amplification master mix have been recorded in the *Chemical Preparation* module of STACS.

Identifiler Direct MM		GlobalFiler Express MM	
	Per Sample (µl)		Per Sample (µl)
IDD Primer Set	12.5	GFE Primer Set	6.0
IDD Master Mix/Enzyme	12.5	GFE Master Mix w/Additive	6.0

- 4.2.3** Within STACS, select plate(s) to be processed and select the appropriate scenario.

- 4.2.3.1** Additionally, for automated processing only:

- Scan the instrument barcode on the Robotic Workstation.
- Ensure the Robot Maintenance Checks have been performed.
- Indicate whether each check passed.

- 4.2.4** If not previously performed during the selection of the plate(s), scan the barcode on each of the 96-Well plate(s). Scan the barcode of each reagent required for the selected scenario. Select "*Process*" and proceed with the Master Mix Addition procedure.

- 4.2.4.1** Additionally, for automated processing only:

- STACS launches the robotic software
- If necessary, enter the appropriate user name and password at the robotic software log-in screen
- Verify that the appropriate Master Mix Addition script has been opened.
- Ensure the required reagents and the selected 96-Well plate(s) have been loaded on to the instrument. On the Tecan EVO, the plates are loaded in the hotels with well A1 in the back right corner.
- Ensure the instrument has been properly flushed and no air bubbles are visible in the tubing or syringes
- Initiate the robotic software process.
- Indicate the number of plates to be processed.
- Start the Master Mix Addition script.

The following Master Mix Addition procedure will be performed manually or by the Robotic Workstation:

- 4.2.5** Add the required reagents to the appropriate wells in the 96-Well plate(s):

4.2.5.1 For Identifiler Direct:

1. Add 25 µL of Identifiler Direct master mix to each allocated well.
 - IDD Ladder wells do not receive master mix.

4.2.5.2 For GlobalFiler Express:

- Add 12 µL of GlobalFiler Express Amplification master mix to each allocated well.
- GFE Ladder wells do not receive master mix.

4.2.6 For automated processing, the Robotic Workstation will heat-seal the 96-Well plate(s) with a foil cover and centrifuge each plate(s) at 1000 RPM for 1 minute. For manual processing, heat-seal the 96-Well plate(s) with a foil cover or seal the plate with an optical cover.

4.2.7 Visually inspect the plate(s) and indicate the result in STACS as successful, failed or aborted. Comments and observations must be entered for plates with process failed results. If the plate(s) were processed on a Robotic Workstation, indicate in STACS whether the bleach process was performed.

4.3 Identifiler® Direct , or GlobalFiler® Express PCR Amplification

4.3.1 If necessary (i.e., manual Master Mix Addition), quick-spin the 96-Well plate(s) in a centrifuge (approximately 30 seconds).

4.3.2 Load each plate into a thermal cycler and place an ABI Optical Compression Pad on each plate and close the thermal cycler(s) lid(s).

4.3.3 Select the appropriate method on the thermal cycler. Available methods are:

Identifiler Direct PCR Amplification Method (on 9700)			
	Temperature (C)	Time	Sample Type and Cycles
Hold	95	11 minutes	
Cycle	94	20 seconds	Buccal (25-28 cycles)
	59	2 minutes	
	72	1 minute	
Hold	60	25 minutes	
Hold	4	Forever (∞)	

GlobalFiler Express PCR Amplification Method (on 9700 or ProFlex)			
	Temperature (C)	Time	Sample Type and Cycles
Hold	95	1 minute	
Cycle	94	3 seconds	Buccal or Blood (26 or 28 cycles)
	60	30 seconds	
Hold	60	8 minutes	
Hold	4	Forever (∞)	

- 4.3.4** Start the thermal cycler.
- **Identifiler Direct:** Ensure the reaction volume is 25 µL and the ramp speed is 9600 on a 9700.
 - **GlobalFiler Express:** Ensure the reaction volume is 15 µL and the ramp speed is MAX on a 9700 or GeneAmp™PCR System 9700 simulation mode on a ProFlex.
- 4.3.5** Ensure the *Thermal Cyclers Bar Code* and the *Plate Bar Code* for each plate to be amplified has been scanned into STACS.
- 4.3.6** Indicate the result of the process in STACS as successful, failed or aborted. Comments and observations must be entered for plates with process failed results.
- 4.3.7** After amplification, store plate(s) refrigerated in the post-amplification laboratory at 4°C ± 3°C.

NOTE: The MicroAmp support base used to transport the 96-Well plate from the pre-amplification laboratory space to the post-amplification laboratory space must be cleaned with a 10% bleach solution and/or with 70% isopropyl alcohol prior to or immediately upon return to pre-amplification laboratory space. This practice will minimize the transport of amplified DNA product from post-amplification laboratory space to pre-amplification space.

5 Sampling

Not applicable.

6 Calculations

Not applicable.

7 Measurement Uncertainty

Not applicable.

8 Limitations

The appropriate processing methods are selected for a plate based on the sample type added to the plate and the amplification kit to be used. Based on internal studies, only the combinations of processes listed below are approved for use.

Sample Type	Amp Kit	MM Addition	Punch Size	Thermal Cycler	Cycles
Buccal (FTA)	Identifiler Direct	Manual or Automated using a Tecan EVO Robotic Workstation	1.2mm	9700	25, 26, 27 or 28
Buccal (FTA or Non-FTA) or Blood (FTA)	GlobalFiler Express	Manual or Automated using a Tecan EVO Robotic Workstation	1.2mm	9700 or ProFlex	26 or 28

9 Safety

9.1 All FDDU samples that contain blood are considered potentially infectious regardless of the perceived status of the source individual or the age of the material. All FDDU personnel who work with such material will follow the “Bloodborne Pathogen Exposure Control Plan” found in the most current version of the *FBI Laboratory Safety Manual*.

9.2 Refer to the “Safe Work Practices and Procedures”, “Bloodborne Pathogen Exposure Control Plan”, “Personal Protective Equipment”, and “Chemical Hygiene Plan” sections of the *FBI Laboratory Safety Manual* for important personal safety information prior to conducting these procedures.

9.3 Refer to the “Hazardous Waste Disposal” section of the *FBI Laboratory Safety Manual* for important information concerning proper disposal of the chemicals used in these procedures as well as the biohazardous wastes generated.

10 References

The procedures described here are derived from a variety of sources. Portions of the protocol come directly from some of the references cited below.

Federal Bureau of Investigation. Quality Assurance Standards for DNA Databasing Laboratories, current version.

FBI Laboratory Quality Assurance Manual.

FBI Laboratory Safety Manual.

DNA Procedures Manual

National Research Council. *The Evaluation of Forensic DNA Evidence*. Washington, DC, National Academy Press (1996).

STACS DNA Inc. *Sample Tracking and Control System (STACS) User's Guide*.

Applied Biosystems. *GeneAmp® PCR System 9700 User's Manual Set*. Applied Biosystems, Foster City, CA. 1997.

Applied Biosystems. *Proflex™ PCR System User Guide*. Life Technologies Corporation | Carlsbad, CA. 2016.

Applied Biosystems. *AmpFlSTR® Identifiler® Direct PCR Amplification Kit Product Insert*, P/N 4408617, Applied Biosystems, Foster City, CA.

Applied Biosystems. *AmpFlSTR® Identifiler® Direct PCR Amplification Kit User Guide*. Applied Biosystems, Foster City, CA.

Applied Biosystems. *GlobalFiler® Express PCR Amplification Kit User's Manual*. Applied Biosystems, Foster City, CA.

Applied Biosystems. *GlobalFiler® Express PCR Amplification Kit User Guide*. Applied Biosystems, Foster City, CA.

Borys, S.; Bowen, K.L; and Fourney, R. PCR Volume Reduction Study Using Bloodstained FTA Collection Cards and Capillary Electrophoresis. 9th International Symposium on Human Identification, Oct. 7-10, 1998.

Budowle, B. STR allele concordance between different primer sets: a brief summary. *Profiles in DNA* (2000) 3: 10-11.

Budowle, B., Masibay, A., Anderson, S.J., Barna, C., Biega, L., Brenneke, S., Brown, B.L., Cramer, J., DeGroot, G.A., Douglas, D., Duceman, B., Eastman, A., Giles, R., Hamill, J., Haase, D.J., Janssen, D.W., Kupferschmid, T.D., Lawton, T., Lemire, C., Llewellyn, B., Moretti, T., Neves, J., Palaski, C., Schueler, S., Sgueglia, J., Sprecher, C., Tomsey, C., and Yet, D. STR primer concordance study. *Forensic Science International* (2001) 124: 47-54.

Del Rio, S.A. Cost-Effectiveness in Sample Processing using the FTA™ Treated Stain Card for High Throughput. 8th International Symposium on Human Identification, Oct. 1997.

DNA Purification Using FTA™ Cards, Extraction of Peripheral Blood DNA for PCR Analysis. Life Technologies Product Information Publication.

Moretti, T.R., Baumstark, A.L., Defenbaugh, D.A., Keys, K.M., Smerick, J.B., and Budowle, B. Validation of short tandem repeat (STRs) for forensic usage: Performance testing of fluorescent

multiplex STR systems and analysis of authentic and simulated forensic samples, *Journal of Forensic Science* (2001) 46: 647-660.

Moretti, TR., Koons, BW, and Budowle, B. Enhancement of PCR amplification yield and specificity using AmpliTaq Gold™ DNA Polymerase, *BioTechniques* (1998) 25: 716-722.

Procedure for the amplification of DNA from blood spotted on FTA™ paper. FTA paper pamphlet. Fitzco, Inc. 1996.

Thacker, C., Phillips, C.P., and Synercombe-Court, D. (2000) Use of FTA cards in small volume PCR reactions. *Progress in Forensic Genetics* 8 (1999). 8, pp. 473-475.

Tindall, KR. and Kunkel, TA. Fidelity of DNA synthesis by the *Thermus aquaticus* DNA polymerase, *Biochemistry* (1988) 27: 6008-6013.

Saiki, RK., Scharf, S., Faloona, F., Mullis, KB., Horn, FT., Erlich, HA. and Arnheim, N. Enzymatic amplification of β -globin sequences and restriction site analysis for diagnosis of sickle cell anemia, *Science* (1985) 230:1350-1354.

Wang et al. Development and Validation of the AmpFlSTR® Identifiler® Direct PCR Amplification Kit: Multiplex Assay for the Direct Amplification of Single-Source Samples, *Journal of Forensic Science* (2011) doi: 10.1111/j.1556-4029.2011.01757.x.

Rev. #	Issue Date	History
8	03/11/16	Complete rewrite of procedure for simplification, elimination of redundancy and to remove software interface instructions. Added FTA extraction steps from 304-7. Moved QA/QC guidance to Appendices.
9	12/09/16	Updated to incorporate implementation of GlobalFiler Express and Non-FTA samples. Simplified Kit QC set-up instructions in Appendix B.
10	09/13/19	Deleted Identifiler, Agilent BioCel, Tecan TE-MO and FTA procedures throughout. Added Proflex ThermalCycler and corresponding parameters. Added SwabSolution Procedures. Appendix A: Removed Ethanol Wash.

Approval

Redacted - Signatures on File

DNA Technical Leader

Date: 09/12/2019

FDDU Chief

Date: 09/12/2019

QA Approval

Quality Manager

Date: 09/12/2019

Appendix A: Instrument Maintenance

Refer to the DNA procedure for equipment calibration and maintenance for additional information and applicable frequency requirements.

1. Robotic Workstations

A. Tecan Instrument Cleaning and General Maintenance Guidance

- **Seal Integrity Check** - Verify that the tubing and syringes (plunger lock screws) are tight and no air bubbles are being introduced into the system.
- **Tip Check** - Verify that the LIHA tips (as appropriate) are tight, free of clogs and not dripping system liquid (RO purified water).
- **ROMA Check** - Verify that the Tecan ROMA(s) and LIHA (as appropriate) are moving and working properly.
- **System liquid level check** - Check and refill the system liquid (i.e., water), if necessary. When refilled, the system liquid should be allowed to de-gas overnight.
- **Clean tips with alcohol** - Clean the outside of the fixed Tecan tips with 70% isopropyl alcohol each workday before use and if they become visually soiled.
- **Bleach deck** - Decontaminate the Tecan workdeck with 10% Bleach each workday before use and if it becomes visually soiled. The barcode scanners should be cleaned with a lint-free cloth only, as needed.
- **System liquid flush** - Flush the instrument thoroughly until no air bubbles are visible in the tubing or syringes.
- **RoboScrub** - Weekly, as used, The Tecan Robotic workstations that utilize the Liquid Handling Arm (LIHA) should be flushed with a mild detergent (e.g., Liquinox) to clean the inside of the Tecan tips, tubing and syringes to maintain the precision and accuracy of liquid handling by the Tecan Robotic Workstation. The RoboScrub solution should remain in the Tecan for a minimum of 15 minutes. Following the elapsed time, flush the Tecan with water until there are no visible air bubbles in the tubing or syringes
- **Bleach process** - Following each run on the Tecan EVO workstations, decontaminate the inside and outside of the Tecan tips with 3% Bleach and flush with water. There is a Bleach Wash script that will perform this process.

B. Performance Verification (aka Artel)

An Artel MVS (Multichannel Verification System) and NIST traceable standards are used to test the accuracy and precision of the liquid handling by a Tecan Robotic Workstation. Refer to the Artel MVS Multichannel Verification System User Guide for operation of the Artel MVS.

For the LIHA (typically configured with eight (8) fixed tips), perform a minimum of 6 repetitions with each tip for each volume.

The accuracy and precision results must be within the tolerance limits set by the FDDU for each volume.

- At times, it may be necessary to modify/optimize the Tecan liquid class parameters (e.g., offset and factor) to ensure that the accuracy and precision of the Tecan continually meets the FDDU tolerance limits.
- The tolerance limits, Artel data, and any modifications to Tecan liquid classes must be maintained in a notebook or an electronically accessible file/folder for reference.

2. Thermal Cycler General Maintenance and Performance Verification

A. 9700

Refer to the *GeneAmp®* PCR System 9700 User's Manual Set 96-Well Sample Block Module User's Manual for instructions on how to perform the following procedures.

- **Cleaning** - Refer to *Cleaning the sample wells* and *Cleaning the sample block cover* sections.
Temperature Verification Test - This test procedure requires the use of a Temperature Verification System (Applied Biosystems) and verifies that the thermal cycler remains within the temperature accuracy specification. Refer to *Running the Calibration Verification Test* section.
Temperature Non-uniformity Test - This test procedure requires the use of a Temperature Verification System and verifies the temperature uniformity of the sample wells in the thermal cycler. Refer to *Running the Temperature Non-uniformity Test* section.
- **Rate Test and Cycle Test** - These procedures verify the integrity of the cooling and heating system of a thermal cycler. Refer to *Running System Performance Diagnostics* section.

B. ProFlex™

Refer to the ProFlex™ PCR System User Guide for instructions on how to perform the following procedures.

- **Cleaning** – Refer to the chapter for maintaining the instrument for instructions on how to clean the sample wells and heated cover.
- **Verify Block Temperature** – These tests are found in the *Block Verification Test* screen and requires the use of a Temperature Verification System for the following test types:
 - **Heated Cover Test** - This test verifies the proper functioning of the heated cover.
 - **Temperature Verification Test** - This procedure verifies that the thermal cycler remains within the temperature accuracy specification.
 - **Temperature Non-Uniformity Test** - This procedure verifies the temperature uniformity of the sample wells in the thermal cycler.

Appendix B: QC of Critical Reagents

Identifiler Direct Kit

Each new lot of Identifiler Direct Kits will be evaluated by testing it jointly with the current lot using the following samples per lot, at minimum:

- 5 - Punches (1.2mm) of a BIS (buccal) control,
- 1 - Combo control,
- 5 - 9947A controls with a template quantity of 4 ng

Amplify the plate at 25 cycles using the Identifiler Direct protocol.

Inject the plate, with at least 2 allelic ladders per lot, 3x at the current CE settings.

Acceptance Criteria:

The allelic ladders from both the new lot and current lot of kits must be used separately to analyze the QC data with GeneMapper ID-X. Criteria for acceptance of kits for use with FDDU samples are:

1. BIS and 9947A samples:
 - a. Acceptable peak morphology and balance at each locus.
 - b. Average peak heights of all alleles are greater than 150 RFU.
 - c. Correct typing results obtained.
 - d. No allelic peaks, other than those attributable to the BIS and 9947A controls, are detected.
2. Each ladder allele is greater than 100 RFU.
3. Acceptable results for the negative controls.

If all the above criteria are met then the new lot is approved for use and the allelic ladders included in both the new lot and current lot of kits may be used for the analysis of samples both within and between lots of Identifiler Direct kits. Enter the results of the QC procedures in STACS using the *Receiving* modules.

GlobalFiler Express Kit

Each new lot of GlobalFiler Express Kits will be evaluated by testing it jointly with the current lot using the following samples per lot, at minimum:

- 5 - punches (1.2mm) of a BIS (buccal) control
- 1 - Combo control
- 4 - 007 controls with a template quantity of 5 ng

Amplify the plate at 26 cycles using the GlobalFiler Express protocol.

Inject the plate, with at least 2 allelic ladders per lot, 3x at the current CE settings.

Acceptance Criteria:

The allelic ladders from both the new lot and current lot of kits must be used separately to analyze the QC data with GeneMapper ID-X. Data will be analyzed using normalization. Criteria for acceptance of kits for use with FDDU samples are:

1. BIS and 007 samples:
 - a. Acceptable peak morphology and balance at each locus.
 - b. Average peak heights of all alleles are greater than 175 RFU.
 - c. Correct typing results obtained.
 - d. No allelic peaks, other than those attributable to the BIS and 007 controls, are detected.
2. Each ladder allele is greater than 175 RFU.
3. Acceptable results for the negative controls.

If all the above criteria are met then the new lot is approved for use and the allelic ladders included in both the new lot and current lot of kits may be used for the analysis of samples both within and between lots of GlobalFiler Express kits. Enter the results of the QC procedures in STACS using the *Receiving* modules.