

Technical Procedures for Processing with Physical Developer

1 Scope

Physical Developer is used by FBI Laboratory Friction Ridge Discipline personnel to develop latent prints on porous, semi-porous, and certain non-porous surfaces.

2 Limitations

All metal items, such as staples and paper clips, must be removed from item(s) prior to Physical Developer processing. Metal tweezers cannot be used during processing.

3 Equipment/Materials/Reagents

Citric Acid

Maleic Acid

Silver Nitrate

Bleach

Synperonic N or Tween 20

Ferrous Ammonium Sulfate

n-Dodecylamine Acetate

Ferric Nitrate

4 Procedures

4.1 Solution Preparation

Personnel will prepare the solutions as follows. Alternative amounts may be prepared, provided the same ratio of chemicals mixed is retained.

4.1.1 Maleic Acid solution

Combine:

- Maleic Acid - 25 g
- Distilled water - 1000 ml

Stir until solid dissolves.

4.1.2 Redox solution

In the order listed below, stir until each solid is dissolved before adding the next solid:

- Distilled water - 1000 ml
- Citric Acid - 20 g
- Ferric Nitrate - 30 g
- Ferrous Ammonium Sulfate - 80 g

4.1.3 Detergent solution

Combine:

- n-Dodecylamine Acetate - 3 g
- Synperonic N or Tween 20- 4 g
- Distilled water - 1000 ml

Stir until all chemicals dissolve.

4.1.4 Silver Nitrate solution

Combine:

- Silver Nitrate - 200 g
- Distilled water - 1000 ml

Stir until solid dissolves.

4.1.5 Physical Developer working solution

While stirring, combine in the order listed:

- Redox solution - 1000 ml
- Detergent solution - 40 ml
- Silver Nitrate solution - 50 ml

Stir for at least three minutes.

4.1.6 Bleach solution (optional)

Combine:

- Bleach - 500 ml
- Distilled water - 500 ml

4.2 Application

Personnel will complete the following steps in order:

1. Immerse item(s) in Maleic Acid solution.
2. Agitate solution, manually or with orbital shaker, for a minimum of 15 minutes.
3. Immerse item(s) in Physical Developer working solution.
4. Agitate solution, manually or with orbital shaker, for 10-15 minutes.
5. Immerse item(s) in first water rinse for at least 1 minute.
6. Rinse item(s) in second water rinse.
7. Dry item(s) in air or by applying heat with an iron, heater, or dryer.

4.2.1 Bleach Solution Rinse (Optional)

The bleach solution may darken latent prints developed with Physical Developer, lighten the background, and remove any previous processing stains that may still be present on the item. The bleach solution is especially effective on paper bags and paper currency. The solution is usually applied after the second water rinse but can be done after the item has been dried.

Personnel will complete the following steps in order:

1. Immerse item(s) in Bleach solution for 10-15 seconds or until desired contrast is achieved.
2. Rinse item(s) in water.
3. Dry item(s) as described above.

4.2.2 For digital capture and photography, see FBI Friction Ridge Discipline Processing Manual Preamble.

4.3 Storage

Maleic Acid, Redox, Detergent, and Bleach solutions may all be stored in any type of laboratory acceptable receptacle.

Silver Nitrate solution must be stored in a dark bottle.

Physical Developer working solution is not stored. It is prepared as needed.

4.4 Shelf Life

Maleic Acid, Redox, and Bleach solution have indefinite shelf lives provided the reagent checks are satisfactory.

Detergent and Silver Nitrate solution each have a shelf life of 1 year provided the reagent checks are satisfactory.

Physical Developer working solution is not retained. It is prepared as needed.

5 Standards and Controls

See FBI Friction Ridge Discipline Processing Manual, Preamble.

6 Safety

See FBI Laboratory Safety Manual for appropriate information.

7 Sampling

Not applicable.

8 Calculations

Not applicable.

9 Measurement Uncertainty

Not applicable.

10 References

Cantu, A. A. "Silver Physical Developers for the Visualization of Latent Prints on Paper." *Forensic Sci Rev.* 13:29.

De la Hunty M, Moret S, Chadwick S, Lennard C, Spindler X, Roux C. "Understanding Physical Developer (PD): Part I – Is PD targeting lipids?" *FSI.* 257:481-487.

De la Hunty M, Moret S, Chadwick S, Lennard C, Spindler X, Roux C. "Understanding Physical Developer (PD): Part II – Is PD targeting eccrine constituents?" *FSI.* 257: 488-495.

FBI Laboratory Safety Manual, Federal Bureau of Investigation, Laboratory Division. Latest Revision.

FBI Friction Ridge Discipline Processing Manual, Preamble, Federal Bureau of Investigation, Laboratory Division. Latest Revision.

Frazier, A. "Optimization of the Physical Developer Process" Honors Intern Project Paper, Federal Bureau of Investigation, Laboratory Division, Quantico VA, 2009.

Houlgrave S, Andress M, Ramotowski, R. "Comparison of Different Physical Developer Working Solutions - Part I: Longevity Studies". *JFI.* 61(6): 621-639.

Houlgrave S, Ramotowski, R. Comparison of Different Physical Developer Working Solutions - Part II: Reliability Studies. *JFI.* Nov/Dec 2011; 61 (6): 640-651.

Phillips, C. E., Cole, D. O., and Jones, G. W. "Physical Developer: A Practical and Productive Latent Print Developer," *JFI.* 40(3):135.

Sauzier G, Frick AA, Lewis, SW. "Investigation into the Performance of Physical Developer Formulations for Visualizing Latent Fingerprints on Paper". *JFI.* 63(1): 70-89.

Trozzi, T. A., Schwartz, R. L., and Hollars, M. L. *Processing Guide for Developing Latent Prints*, FBI Laboratory, Washington DC, 2001.

Rev. #	Issue Date	History
2	10/02/17	Updated for Biometrics Analysis Unit. Section 5.3, modified “accepted” to “acceptable”. Abbreviations addressed.
3	07/15/21	Replace Latent Print Units with Friction Ridge Discipline. Minor wording changes. Streamline equipment list. Re-organization and re-numbering of sections. Section 4.4 broken into Section 4.1.1, Section 4.1.2, Section 4.1.3, Section 4.1.4, Section 4.1.5, and Section 4.1.6. Parts of Section 4.2 broken into Section 4.2.1 and Section 4.2.2. Section 4.1, added ratio allowance. Section 4.1.6, moved note. Section 4.2, removed note. Section 5, added Preamble

Approval

Redact - Signatures on File

Friction Ridge Discipline
Technical Leader

Date: 07/14/2021

Latent Print Operations
Unit Chief

Date: 07/14/2021

Latent Print Support Unit
Chief

Date: 07/14/2021

Scientific and Biometrics
Analysis Unit Chief

Date: 07/14/2021